

Thermo Scientific Dionex AutoTrace 280 Solid-Phase Extraction Instrument

The Thermo Scientific™ Dionex™ AutoTrace™ 280 Solid-Phase Extraction instrument handles large volume Solid-Phase Extraction (SPE) of liquid samples from 20 mL to 4 liters. The instrument can be configured with 1, 3, or 6 mL SPE cartridges or 47 mm solid-phase disks. The instrument automatically loads liquid samples on the SPE material then automatically extracts the analytes of interests using up to five different solvents.



System Features

- Automatically conditions, rinses, and elutes SPE cartridges with a choice of five solvents
- Simultaneously loads up to six samples
- Loads sample volumes from 20 mL to 4 liters
- Offers choice of six different collection vial racks
- Segregates aqueous and solvent waste using separate waste lines
- Processes ground, surface, and wastewater samples
- Easy-to-use SPE software
- Supports SPE cartridges
- Documented EPA methods for drinking, surface, and waste water
- Integrated venting of solvent vapors

Applications

The Dionex AutoTrace 280 SPE instrument is perfect for laboratories looking to automate large volume water extractions. The Dionex AutoTrace 280 instrument applications include many of the U.S. EPA methods and a variety of other SPE sample preparation cleanup methods.

Methods

- U.S. EPA Method 508 & 608
- U.S. EPA Method 525
- U.S. Method 539
- U.S. EPA Method 549.2
- U.S. EPA Method 625
- Extraction of organochlorine pesticides
- Extraction of endocrine disruptors
- Extraction of dioxins and furans
- Extraction of triclosan

Table 1. Shows a pesticide recovery study comparing the AutoTrace SPE technique and a vacuum manifold technique. The improvements in recovery and reproducibility are attributed to the microprocessor control of all the liquid flow rates – both sample and SPE solvents.

Pesticide Recovery Study N = 6	Dionex AutoTrace 280 SPE		Vacuum Manifold SPE	
	Recovery%	RSD	Recovery%	RSD
Compound				
Atrazine	88	1.8	54	12.2
Propazine	91	1.5	80	7.3
Alachor	99	3.4	96	4.1
Metachor	99	4.3	96	2.9

- Extraction of nitrosamines

SPECIFICATIONS

Gas Regulator and Gas Gauge Range	Output: 0–30 psi (0–1.4 bars), Input: 100 psi (6.9 bars) maximum
Net Weight	150 lbs. (68.1 kg)
Dimensions (h × w × d)	23 × 25 × 27 in (57 × 63.5 × 69 cm) Sample Rack: 8 × 16.5 × 13 in)
500 & 600 Series EPA Methods	24 methods stored in the Dionex AutoTrace 280 SPE instrument
Electrical	Voltage: 100, 120, 220, or 240 V ± 10% Frequency: 47–63 Hz Power: 150 Volt AC
Liquid Management	Air Syringe: One 10 mL air syringe LH Syringe: One 10 mL liquid handling syringe 12-port Valve: Rotary, sliding Rulon® seal Valves: 3-way, Teflon® Nozzles: Stainless steel Tubing: TFE, 1/16" ID, (6 total supplied)
Sample Pumps	Displacement: Positive Accuracy: ± 2.5% Tube Fitting: Kynar® Piston and Liner: Ceramic Non Use: Not for use with acetic acid and acetone
Solid Phase Extraction Configurations	1 mL Syringe: Compatible cartridge 3 mL Syringe: Compatible cartridge 6 mL Syringe: Compatible cartridge 47 mm SPE disk

Ordering Information

In the U.S., call (800) 346-6390 or contact the Thermo Fisher Scientific Regional Office nearest you. Outside the U.S., order through your local Thermo Fisher Scientific office or distributor. Refer to the following part numbers.

Description	Part Number
Dionex AutoTrace 280 Automated Large Volume SPE for Disks	071386
Dionex AutoTrace 280 Automated Large Volume SPE for 6 mL Cartridges	071385
Dionex AutoTrace 280 Automated Large Volume SPE for 3 mL Cartridges	072605
Dionex AutoTrace 280 Automated Large Volume SPE for 1 mL Cartridges	072604
Dionex AutoTrace 280 Automated Large Volume SPE for 6 mL Glass Cartridges	072606

www.thermoscientific.com/samplepreparation

©2014 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. Rulon is a registered trademark of Saint Gobain. Teflon is a registered trademark of E.I. du Pont de Nemours and Company. Kynar is a registered trademark of Arkema Inc. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Brazil +55 11 3731 5140
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 6 6885 1213
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
Taiwan +886 2 8751 6655
UK/Ireland +44 1442 233555
USA +1 800 532 4752

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand