

**Thermo Scientific Reagents, Solvents and Accessories** 

- Gas Chromatography
- High Performance Liquid Chromatography
- Mass Spectrometry





# Thermo Scientific Products for GC and HPLC

Derivatization Reagents for Specific Functional Groups	s 1-2
Derivatization Reagents for Drug-of-Abuse	3-4
Troubleshooting	5
Gas Chromatography Reagents 6-25	5
Silylation Reagents	7-14
Acylation Reagents	15-18
Alkylation Reagents	19-21
Siliconizing Fluids	23-24
Books	25
HPLC Ion Pair Reagents 26-31	
HPLC Ion Pair Reagents	26-28
Derivatization and Visualization Reagents for HPLC	29-31

# Amino Acid Analysis 33-43

Introduction to Amino Acid Analysis	33-34
Fluorometric Detection Reagents for Amino Acids	35
Ninhydrin Detection Reagents for Amino Acids	36-37
Sample Preparation and Hydrolysis	38
Hydrolysis Reagents	39
High-Purity Pre-Column Derivatization Reagents	40-41
Solvents for Amino Acid Analysis	42
Peptide Retention Standards	43
Sample Handling 44-51	
Heating/Stirring/Evaporation	44-47
Vials and Closures	48-51
Alphabetical Index	52
Numerical Index	53

inctional Group	Procedure	Reagent	Derivative	Notes
	Silylation	BSA	TMS Amides	Difficult to form due to steric hindrance
	,	BSTFA	TMS Amides	
Amides		BSTFA+TMCS	TMS Amides	TMCS used as a catalyst
		MSTFA	TMS Amides	Reaction byproducts more volatile
0		MSTFA+TMCS	TMS Amides	,,,
II -C-NH <sub>2</sub>		Tri-Sil Reagents	TMS Amides	
-		MTBSTFA	TBDMS Amides	Difficult to form; very stable
Primary		MTBSTFA+TBDMCS	TBDMS Amides	TBDMCS aids derivatization
	Al - +i			I DDIVICS alus derivadization
0	Acylation	MBTFA	Trifluoroacetamides	
<u> </u>		TFAA	Trifluoroacetamides	
-C-NHR		PFAA	Pentafluoropropionamides	Good for ECD detection
Secondary		HFBI	Heptafluorobutyamides	
	Alkylation	MethElute Reagent (TMPAH)	Methyl Amides	On-column derivatization especially for drugs
	Silylation	BSA	TMS	ioi urugs
		BSTFA	TMS	
Amines		BSTFA+TMCS	TMS	TMCS aids derivatization
		MSTFA	TMS	
H		MSTFA+TMCS	TMS	TMCS aids derivatization
-C-NH <sub>2</sub>		Tri-Sil® Reagents	TMS	
-0-1411 <sub>2</sub>	Silylation	MTBSTFA	TBDMS	Difficult to form, but more stable
Н	Silyiduoii	MTBSTFA+TBDMCS	TBDMS	TBDMCS aids derivatization
Primary	Aculation		Trifluoroacetamides	
	Acylation	MBTFA	Trifluoroacetamides	Good for trace analysis with ECD
H		TFAA		Good for trace analysis with ECD
–Ċ–NHR		TFAI	Trifluoroacetamides	Good for trace analysis with ECD
ı		PFAA	Pentafluoropropionamides	
Н		PFPI	Pentafluoropropionamides	
Secondary		HFAA	Heptafluorobutyamides	
		HFBI	Heptafluorobutyamides	
	Alkylation	MethElute Reagent (TMPAH)	Methyl Amides	On-column derivatization for specific dru
	Silylation	BSA	TMS	Most often used derivatives
	•	BSTFA	TMS	Good thermal stability
		BSTFA+TMCS	TMS	Poor hydrolytic stability
		HMDS	TMS	Weak donor usually used with TMCS
		MSTFA	TMS	
		MSTFA+TMCS	TMS	
		TMCS	TMS	Weak donor usually used with HMDS;
Hydroxyl-OH				can be used with salts
R-OH		TMSI	TMS	Can be used with syrups
Alcohols		Tri-Sil Reagents	TMS	
		MTBSTFA	TBDMS	More stable than TMS, good MS
<b>⟨○⟩</b> — он		MTBSTFA+TBDMCS	TBDMS	fragmentation patterns TBDMCS aids derivatization
	Acylation	MBTFA	Trifluoroacetates	Good for trace analysis with EDC
Phenols	Acylation	TFAA	Trifluoroacetates	Good for trace analysis with EDC
		TFAI	Trifluoroacetates	Good for trace analysis with EDC
		PFAA		Good for trace analysis with EDC
			Pentafluoropropionates	
		HFBI	Heptafluorobutrates	Good for trace analysis with EDC
		HFAA	Heptafluorobutrates	Good for trace analysis with EDC
	Alkylation	PFBBr	Pentafluorobenzyl Ethers	With alkoxides only
	Silylation	BSTFA	TMS	
Amino		TMSI	TMS	
Acids	Acylation	HFBI	Heptafluorobutrates	
н н он	Alkylation	Methylate Reagent	Methyl Esters	
п ;; оп				

unctional Group	Procedure	Reagent	Derivative	Notes
Catecholamines	Silylation	TMSI	TMS	
ОН	Acylation	HFAA PFAA TFAA HFBI	Heptafluorobutrates Pentafluoropropionates Trifluoroacetates Heptafluorobutrates	
	Silylation	MSTFA	TMS	
Carbohydrates	Silylation	TMSI Tri-Sil Reagents	TMS TMS	Can be used with some syrups
(CH <sub>2</sub> OH) <sub>n</sub>	Acylation	MBTFA TFAI	Trifluoroacetates Trifluoroacetates	Volatile derivatives of mono-, di- and trisaccharides
	Silylation	BSA	TMS	Easily formed, generally not stable, analyze quickly
Carboxyl		BSTFA BSTFA+TMCS MSTFA TMCS TMSI Tri-Sil Reagents	TMS TMS TMS TMS TMS TMS TMS TMS	Can be used with some salts
0     -C-0H		MTBSTFA	TBDMS	More stable than TMS derivatives
-с-ип		MTBSTFA+TBDMCS	TBDMS	TBDMCS aids derivatization
	Alkylation	PFBBr BF <sub>3</sub> -Methanol Methylate Reagent (DMFDMA) MethElute Reagent (TMPAH) PFAA+Pentafluoropropanol	Pentafluorobenzyl Esters Methyl Esters Methyl Esters Methyl Esters Methyl Esters Pentafluoropropyl Ester	Used in EC detection & UV, MS  Best for large samples of fatty acids Fatty acids and amino acids On-column derivatization  Drug analysis
Inorganic Anions	Silylation	BSTFA MTBSTFA	TMS TBDMS	,
Nitrosamines	Acylation	HFAA	Heptafluorobutrates	
R <sup>1</sup> N R <sup>2</sup>   N N N N N N N N N N N N N N N N N N				
Cultonomidee	Silylation	BSTFA	TMS	
Sulfonamides  O O R <sup>1</sup> N R <sup>3</sup>	Acylation	HFAA PFAA TFAA	Heptafluorobutrates Pentafluoropropionates Trifluoroacetates	

Amphetamines Amphetamines	BSTFA HFAA	1
Amphetamines	HFAA/PFAA	2-5 6
Amphetamines Amphetamines Methamphetamine	MSTFA with TMCS TFAA TFAA	7 7,8 9,10
·	BSTFA MethElute Reagent (TMPAH) Methylate Reagent (DMFDMA) PFBBr	1 11-13 14,15 16
Benzoylecgonine	BSTFA/Butyl Iodine/TMPAH BSTFA MTBSTFA PFAA/PFPOH	17 1,18 19 9,20
	BSA BSTFA MSTFA TFAI	21 22 21 23
THC metabolites	BSA BSTFA/BSTFA+1% TMCS BSTFA/TMCS/TMSI MSTFA MSTFA/MSTFA+1% TMCS MTBSTFA PFBBr PFAA/HFIOH PFAA/PFPOH TFAA & BF <sub>3</sub> /MeOH MethElute Reagent (TMPAH)	24 24-27 24 9 27 28 29 30 31 32 9
Morphine	BSTFA+1% TMCS MBTFA PFAA TFAA	33 34 35 36
Morphine/Codeine	BSTFA+1% TMCS BSTFA/TFAI HFBA MBTFA PFAA PFAA/HFAA PFAA/PFPOH TFAA	1,37 38,39 40 38 38 38,41 37 9 42 43
PPC/PCHP/PCP	BSTFA+1% TMCS HFAA	44 45
	Methamphetamine  Benzoylecgonine  THC metabolites  Morphine  Morphine/Codeine	Methamphetamine  BSTFA  MethElute Reagent (TMPAH) Methylate Reagent (DMFDMA) PFBBr   Benzoylecgonine  BSTFA/Butyl lodine/TMPAH BSTFA MTBSTFA MTBSTFA PFAA/PFPOH  THC metabolites  BSA BSTFA/BSTFA+1% TMCS BSTFA/TMCS/TMSI MSTFA MSTFA MSTFA PFBBr  PFAA/HFIOH PFAA/PFPOH  TFAA & BF/MeOH MethElute Reagent (TMPAH) TMSI MOrphine  BSTFA+1% TMCS MBTFA PFBBR MBTFA PFBBR MBTFA PFBBR PFAA/HFIOH PFAA/PFPOH TFAA & BF/MeOH MethElute Reagent (TMPAH) TMSI MOrphine  BSTFA+1% TMCS MBTFA PFAA HFBA MBTFA PFAA PFAA PFAA PFAA PFAA PFAA/HFAA PFAA/HFAB PFAA/

See references on following page.
† Reagent names correspond to product names as listed in this catalog, except PFPOH (pentafluoropropanol).

HFIOH (heptafluoro-isopropanol) is not offered by Thermo Fisher Scientific. PFAA (PentaFluoropropionic Acid Anhydride) and HFAA (HeptaFluorobutyric Acid Anhydride) are sometimes incorrectly referred to as PFPA and HFBA (respectively), which are the appropriate abbreviations for the free acid.

#### References

- 1. Dutt, M.C. (1982). J. Chromatogr. 248, 115-124.
- 2. Wu, A.H.B., et al. (1992). J. Anal. Toxicol. 16, 137-141.
- 3. Thurman, E.M., et al. (1992). J. Anal. Toxicol. 16, 19-27.
- 4. Reimer, M.L., et al. (1993). Biol. Mass Spectrom. 22, 235-242.
- 5. Yamamoto, T., et al. (1989). J. Anal. Toxicol. 13, 117-119.
- 6. Wu, A.H.B., et al. (1992). Biol. Mass Spectrom. 21, 278-284
- 7. Rood, H.D. and Knitter, J.A. (1991). Capillary Chromatography. W.G. Jennings, ed., pp. 115-130.
- 8. DePace, A., et al. (1990). J. Forensic Sci. 35(6), 1431-1435.
- 9. Mulé, S.J. and Casella, G.A. (1998). J. Anal. Toxicol. 12, 102-107.
- 10. Suzuki, S., et al. (1983). J. Chromatogr. 267, 381-387.
- 11. Kananen, G., et al. (1972). J. Chromatogr. Sci. 10, 283-287.
- 12. Skinner, R., et al. (1973). Anal. Chem. 45(3), 574-576.
- 13. Mulé, S.J. and Casella, G.A. (1989). J. Anal. Toxicol. 13, 13-16.
- 14. Christophersen, A.S. and Rasmussen, K.E. (1980). J. Chromatogr. 192, 363-374.
- 15. Venturella, V.S., et al. (1973). J. Pharm. Sci. 62(4), 662-668. 16. Walle, T. (1975). J. Chromatogr. 114, 345-350. 17. Ellerbe, P., et al. (1992). J. Anal. Toxicol. 16, 158-162.

- 18. Graffeo, A.P., et al. (1976). J. Chromatogr. 126, 712-717.
- 19. Harkey, M.R., et al. (1991). J. Anal. Toxicol. 15, 260-265.
- 20. Bodor, G., et al. (1990). Clin. Chem. 36, 742-747.
- 21. Jane, I. and Wheals, B.B. (1973). J. Chromatogr. 84, 181-186.
- 22. Nelson, C.C. and Foltz, R.L. (1992). Anal. Chem. 64, 1578-1585.
- 23. Lim, H.K., et al. (1988). Anal. Chem. 60, 1420-1425.
- 24. Harvey, D.J. (1981). Biomed. Mass. Spec. 8(8). 25. Parry, R.C., et al. (1990). J. Anal Toxicol. 14, 39-43.

- 26. Clatworthy, A.J. (1990). Forensic Sci. Int. 46, 219-230.
- 27. Baker, T.S., et al. (1984). J. Anal Toxicol. 8, 255-259.
- 28. Clouette, R., et al. (1993). J. Anal Toxicol. 17, 1-4.
- 29. Rosenfeld, J.M., et al. (1986). Anal. Chem. 58, 716-721.
- 30. McBurney, L.J., et al. (1986). J. Anal Toxicol. 10, 56-64.
- 31. Karlsson, L., et al. (1983). J. Anal Toxicol. 7, 198-202. 32. Foltz, R.L. (1984). Advances in Analytical Toxicology, I. R.C. Baselt, ed. Foster City, CA: Biomedical Publications.
- 33. Wilkinson, G.R., et al. (1969). Biochem. Pharmacol. 18, 1435-1439.
- 34. Weitz, C.J., et al. (1986). Proc. Nat. Acad. Sci. 83, 9784-9788.
- 35. Fehn, J. and Megges, G. (1985). J. Anal. Toxicol. 9, 134-138.
- 36. Wasels, R., et al. (1989). J. Chromatogr. 489, 411-418.
- 37. Christopherson, A.S., et al. (1987). J. Chromatogr. 422, 117-124.
- 38. Chen, B.H., et al. (1990). J. Anal. Toxicol. 14, 12-17.
- 39. Bowie, L. and Kirkpatrick, P.B. (1989). J. Anal. Toxicol. 13, 326-329.
- 40. Nelson, C.C. and Foltz, R.L. (1992). Anal. Chem. 64, 1578-1585.
- 41. Grinstead, G.F. (1991). J. Anal. Toxicol. 15, 293-298.
- 42. Fuller, D.C. and Anderson, W.H. (1992). J. Anal. Toxicol. 16, 315-318.
- 43. Lee, H.M. and Lee, C.W. (1991). J. Anal. Toxicol. 15, 182-187.
- 44. Woodworth, J.R., et al. (1984). J. Anal. Toxicol. 8, 2-6.
- 45. Cone, E., et al. (1981). J. Chromatogr. 223, 331-339.

# **Drugs-of-Abuse Derivatization Applications**

#### **Cocaine Metabolites**

Cardenas, S., et al. (1996). Rapid Commun. Mass Spectrom. 10, 631-636.

Crouch, D.J., et al. (1995). J. Anal. Toxicol. 19, 352-358.

Smimow, D., et al. (1995). Presented at the Society of Forensic Toxicolologists, Baltimore, MD, 72.

Okeke, C.C., et al. (1994). Chromatographia 38, 52-56.

Cone, E.J., et al. (1994). Clin. Chem. 40, 1299-1305.

Peterson, K.L., et al. (1994). Presented at the 46th Annual Meeting of the American

Academy of Forensic Sciences, San Antonio, TX, 197. Moore, J.M., et al. (1993). J. Forensic. Sci. 38, 1305-1325.

Aderjan, R.E., et al. (1993). J. Anal. Toxicol. 17, 51-55.

Taylor, R.W., et al. (1991). Adv. Lab Autom. Rob. 7, 567-582. Ortuno, J., et al. (1990). J. Pharmaceut. Biomed. Anal. 8, 911-914.

Verebey, K., et al. (1989). J. Forensic Sci. 34, 46-52.

Vasiliades, J. (1989). J. Anal. Toxicol. 13, 127.

Isenschmid, D.S., et al. (1988). J. Anal. Toxicol. 12, 242-245.

Mule, S.J., et al. (1988). J. Anal. Toxicol. 12, 153-155.

#### **Derivatization of Cannabinoids**

Prest, H. Advantages of positive chemical ionization mass spectroscopy. Application Note. Hewlett Packard Co.

Lisi, A.M., et al. (1993). J. Chromatogr. 617, 265-270.

Wu, A.H.B., et al. (1993). J. Anal. Toxicol. 17, 215-217.

Clouette, R., et al. (1993). J. Anal. Toxicol. 17, 1-4.

Fysh, R.R. (1989). Presented at the 26th International Meeting of the International

Association of Forensic Toxicologists, Glasgow, Scotland, 67.

Rosenfeld, J.M., et al. (1989). Anal. Chem. 61, 925-928.

Baker, T.S., et al. (1984). J. Anal. Toxicol. 8, 255-259.

Karlsson, L., et al. (1983). J. Anal. Toxicol. 7, 198-202. Borys, H.K., et al. (1981). J. Chromatogr. 205, 303-323.

#### **Derivatization of Amphetamines**

Kuroda, N., et al. (1998). J. Chromatogr. A. 798, 325-334. Dasgupta, A., et al. (1998). Amer. J. Clin. Pathol. 109, 527-532.

Cheung, S., et al. (1997). J. Chromatogr. B. 690, 77-87.

Sadeghipour, F., et al. (1997). J. Chromatogr. A. 761, 71-78.

Hara, K., et al. (1997). J. Anal. Toxicol. 21, 54-58.

Shin, H.-S., et al. (1996). Anal. Chem. 68, 3015-3020.

Meatherall, R.C. (1994). Presented at the 1994 TIAFT-SOFT Program, Tampa, FL, 66.

Thompson, W.C., et al. (1994). Clin. Chem. 40, 1703-1706.

Sievert, H.J.P. (1994). Chirality 6, 295-301.

Melgar, R., et al. (1993). J. Anal. Toxicol. 17, 399-402.

Jones. J.B., et al. (1993). J. Anal. Toxicol. 17, 447.

Gjerde, H., et al. (1993). J. Anal. Toxicol. 17, 65-68.

Cody, J.T., et al. (1992). J. Chromatogr. 580, 77-95.

Hughes, R.O., et al. (1991). J. Anal. Toxicol. **15**, 256-259. Gan, B.K., et al. (1991). J. Forensic Sci. **36**, 1331-1341.

Lillsunde, P., et al. (1991). Forensic Sci. Int. 49, 205-213.

Rasmussen, S., et al. (1989). J. Anal. Toxicol. 13, 263-267. Czarny. R.J., et al. (1989). J. Anal. Toxicol. 13, 257-262.

Hornbeck, C.L., et al. (1989). J. Anal. Toxicol. 13, 144-149.

#### **Derivatization of Opiates**

Brendler, J.P., et al. (1998). Presented at the AAFS, San Francisco, CA.

Balikova, M., et al. (1998). Forensic Sci. Int. 94, 201-209.

Emara, S. (1998). Biomed. Chromatogr. 12, 15-20.

Hyotylainen, T., et al. (1997). J. Chromatogr. A. 771, 360-365.

Guillot, J.G., et al. (1997). J. Anal. Toxicol. 21, 127-133. Dietzen, D.J., et al. (1995). J. Anal. Toxicol. 19, 299-303.

Hendrickson, T.L., et al. (1994). J. Chromatogr. B. 653, 147-154.

Nakahara, Y., et al. (1992). Arch. Toxicol. 66, 669-674.

Chari, G., et al. (1991). J. Chromatogr. 571, 263-270.

Lee, H.-M., et al. (1991). J. Anal. Toxicol. 15, 182-187. Chen, B.H., et al. (1990). J. Anal. Toxicol. 14, 12-17.

McLean, C.F., et al. (1990). J. Pharm. Pharmacol. 42, 669-671.

Wasels, R., et al. (1989). J. Chromatogr. 489, 411-418. Battah, A.-K., et al. (1989). Presented at the 26th Annual International Meeting

of the International Association of Forensic Toxicologists, Glasglow, Scotland, 16.

Derks, H.J.G.M., et al. (1986). J. Chromatogr. 370, 173-178. Moore, J.M., et al. (1984). Anal. Chem. 56, 642-646.

Derivatization Problem	Possible Cause	Recommended Solution
Low Yield	Carrier, air, detector (FID) hydrogen or make-up gas flow set incorrectly	Measure flows using a Thermo Scientific GFM Pro Gas Flow Meter and set accordingly using instrument manufacturer's recommendations
	Reagent deteriorated	Store reagent properly to prevent oxygen/water contamination, temperature damage (refer to product specification sheet)
	Rate of reaction too slow	Re-evaluate reagent concentration, time, temperature and consider heating the reaction mix (consider the thermal stability of the analytes and reagents)
	Water in reaction mix	Remove water by adding sodium sulfate to sample. Store reagent properly to prevent oxygen/water contamination
	Improper handling technique: (e.g. Low boiling components could be lost during sample concentration); sample too dilute; wrong solvent	Re-evaluate technique, if possible eliminate steps in which analyte could be adsorbed or otherwise lost (unnecessary transfers etc.)
	Wrong reagent	Re-evaluate reagent selection and select more appropriate reagent
	Impurities in solvent, starting material, catalysts, or extract interfering with derivatization (e.g. Plasticizers from vial, inorganics used in sample synthesis, preservatives or antioxidants in solvents)	Use only highest purity material at all steps in the sample preparation process
	Reagent: sample ratio too low	Use more reagent for same amount of sample
	Sample adsorbed to glassware	Deactivate glassware, inlet sleeve and column by silanizing
No sample separation after adding reagent and heating	Septum in reaction vial not sealed	Prepare a new sample and derivatize. Be sure that the vial is sealed
Detector response low	Sample components absorbed by inlet liner or column	Inject standard on column known to be performing well. If results are good, remove inlet liner and check cleanliness. Use new, deactivated liner or replace glass wool and packing. Rinse bonded phase column or remove a few cm from inlet end of non-bonded column. If performance is not restored, replace column
	Low yield of derivative - reaction did not go to completion	Add more reagent, increase temperature or heating time or add catalyst. Water may be present; add sodium sulfate to sample
	Detector (FID) dirty	Clean FID as per instrument manual
Extra peak(s)	Derivative reacting with solvent	Use a solvent that does not have an active hydrogen, alcohol or enolizable ketone group (e.g. Hexane, toluene etc. )
	Impurities from sample, solvent, reagents, sample vial, other labware	Inject solvent and reagent blanks, solvent rinse from unused vial etc. Isolate source of impurities
	Reagents interacting with column	Verify that reagent is compatible with analytical column
	Derivative undergoing hydrolysis	Remove water by adding sodium sulfate to sample. Store reagent properly to prevent oxygen/water contamination
Missing peaks or	Wrong reagent	Re-evaluate reagent selection
solvent peak only	Reagent deteriorated	Store reagent properly to prevent oxygen/water contamination, temperature damage (refer to product specification sheet for storage conditions and reagent lifetime)
	Rate of reaction too slow	Re-evaluate reagent concentration, time, temperature and consider heating the reaction mix (consider the thermal stability of the analytes and reagents)
	Impurities in solvent, starting material, catalysts, or extract interfering with derivatization (e.g. Plasticizers from vial, inorganics used in sample synthesis, preservatives or antioxidants in solvents)	Use only highest purity material at all steps in the sample preparation process
	Sample adsorbed to glassware	Deactivate glassware, inlet sleeve and column by silanizing
	Reagent: sample ratio too low	Use more reagent for same amount of sample
	Water in reaction mix	Remove water by adding sodium sulfate to sample. Store reagent properly to prevent oxygen/water contamination



#### Derivatization

The chemical literature contains an abundance of data on derivatization, most of which is relevant to particular compounds, classes of compounds and derivatization reagents. Two books are recognized as standards in the field of analytical derivatization. The first book, Handbook of Analytical Derivatization Reactions by Daniel R. Knapp<sup>1</sup>, provides a general collection of analytical derivatization methods for chromatography and mass spectroscopy (MS) that involves formation of covalent derivatives prior to analysis. The second book, Silylation of Organic Compounds by Alan F. Pierce,2" was a significant factor in the transfer of silvlation reactions from the relatively esoteric field of organosilicon chemistry to the status of perhaps the most widely practiced of derivatization methods."3

#### What is GC Derivatization?

Derivatization is the process of chemically modifying a compound to produce a new compound which has properties that are suitable for analysis using a GC.

# Why Do we Derivatize?

Compounds or compound mixtures are derivatized before analysis for the following reasons:

- 1. To make a compound that otherwise could not be analyzed by a particular method suitable for analysis4
- 2. To improve the analytical efficiency of the compound<sup>5,6</sup>
- 3. To improve the detectability of the compound7

#### Suitability

Often compounds cannot be analyzed because they are not in a form that is suitable for the particular analytical technique. Examples include nonvolatile compounds for GC analysis, 8,9,10 insoluble compounds for HPLC analysis and materials that are not stable using the conditions of the technique. 11 The derivatization procedure modifies the chemical structure of the compounds, allowing analysis by a desired technique.12

#### Efficiency

Direct analysis can be difficult when compounds interact with each other or with the column. These interactions can lead to poor peak resolution and/or asymmetrical peaks that make proper peak integration difficult or impractical. This interference can be reduced with conversion to derivatized products. 13,14 Compounds that exhibit co-elution can often be separated by using the appropriate derivatization methods.

#### Detectability

As demand increases for the analysis of increasingly smaller amounts of materials, it becomes important to extend the detectability range of the materials in question. This increased sensitivity can be accomplished by improved detector design that is directed toward specific atoms or functional groups.

Another popular approach to increase detectability is the use of derivatization. Enhanced detectability can be achieved by increasing the bulk of the compound, or by introducing atoms or functional groups that strongly interact with the detector. 16,17 This technique is performed in gas chromatographic applications, with the addition of halogen atoms for electron capture detectors, 18,19 and with the formation of TMS derivatives to produce readily identifiable fragmentation patterns and mass ions.20

#### Types of Derivatization

Compounds containing functional groups with active hydrogens (-COOH, -OH, -NH and -SH) are usually derivatized prior to analysis by gas chromatography. These functional groups have a tendency to form intermolecular hydrogen bonds that affect the volatility, their tendency to interact deleteriously with column packing materials and their thermal stability.

The ideal derivatization procedure will:

- 1. Accomplish the desired modification.
- 2. Proceed quantitatively, or at least reproducibly.
- 3. Produce products that are readily distinguishable and separable from the starting materials.
- 4. Proceed rapidly with simple and straight-forward laboratory techniques that will be both selective and applicable to a number of similar compounds.
- 5. Involve reagents and reactions that present no unusual hazards.

#### Main Types of Derivatization

- Silylation
- Acylation
- Alkylation

#### Silylation and Silylation Reagents

#### Only Thermo Scientific Reagents offer the combination of variety, quality and reliability.

Silvl derivatives are the most widely used derivatives for gas chromatographic applications. Usually they are formed by the replacement of the active hydrogens from acids, alcohols, thiols, amines, amides and enolizable ketones and aldehydes with the trimethylsilyl group. A variety of reagents is available for the introduction of the trimethylsilyl group. These reagents differ in their reactivity, selectivity and side reactions and the character of the reaction products from the silvlation reagent itself. Considerable literature is available to assist you in the selection of the most suitable silvlation reagent for your particular compounds or systems.1,2

Silvlation reagents and trimethylsilyl derivatives are hydrolytically unstable and must be protected from moisture. However, the rate of hydrolysis for various reagents and derivatives is different, and sometimes it is possible to prepare derivatives in the presence of small amounts of moisture,<sup>21</sup> or to isolate and purify derivatives by extraction in an organic solvent, followed by washing with aqueous solutions.22 Reagents that introduce a t-butyldimethylsilyl group instead of the trimethylsilyl group were developed for greater hydrolytic stability.23 These derivatives provide improved stability against hydrolysis and provide distinctive fragmentation patterns, making them useful in GC/MS applications.<sup>24</sup>

Most trimethylsilyl and t-butyldimethylsilyl derivatives offer excellent thermal stability and are suitable for a wide range of injector and column conditions. However, as the silylation reagents will derivatize nearly all active hydrogens, it is important that they are not injected onto any column in which the stationary phase contains these functional groups. Examples of packings that are not compatible with silylating reagents are polyethylene glycols (TG-WaxMS, TR-WAX or TR-WaxMS) and free fatty acid phases (TR-FFAP).

#### References

- 1. Knapp D.R. (1979). Handbook of Analytical Derivatization Reactions, John Wiley & Sons: New York
- 2. Pierce, A.E. (1968). Silylation of Organic Compounds, Pierce Chemical: Rockford, IL
- 3. Pierce, A.E. (1968). Silylation of Organic Compounds, Pierce Chemical: Rockford, IL. p. 2.
- 4. Sweeley, C.C., et al. (1963). Gas-liquid chromatography of trimethylsilyl derivatives of sugars and related substances. J. Am. Chem. Soc. 85,
- 5. Khalifa, S. and Mumma. R.O. (1972). J. Agric. Food Chem. 20, 632.
- 6. Sinsheimer, J.E. and Smith, R.V. (1967). Methods for the qualitative and quantitative analysis of some hydroxystilbenes. J. Pharm. Sci. 56, 1280.
- 7. Poole. C.F. (1976). Chem. Ind. (London), 479.
- 8. Brittain, G.D. and Schewe, L. (1971). In Recent Advances in Gas Chromatography, Domsky, I.I. and Perry, J.A., Eds., Marcel Decker: New York NY
- 9. Sakauchi, N. and Horning, E.C. (1971). Anal. Lett. 4, 41-42.
- 10. Sullivan, J.E. and Schewe, L.R. (1977), Preparation and gas chromatography of highly volatile trifluoroacetylated carbohydrates using N-Methylbis (trifluoroacetamide) J. Chromatogr. Sci. **15**, 196-197.
- 11. Hallgren, B. and Larsson, S. (1962). J. Lipid Res. 3, 31.
- 12. Samuelsson, K. and Samuelsson, B. (1969). Gas-liquid chromatography-mass spectroscopy of cerebrosides as trimethylsilyl ether derivatives. Biochem. Biophys. Res. Commun, 37(1), 15-21.
- 13. Langer, M., et al. (1960). Chem. Ind. (London), 378.
- 14. Neeman, M., et al. (1959). Tetrahedron, 6, 36.
- 15. Marfey, P. (1984). Determination of D-amino acids. II. Use of a bifunctional reagent, 1,5-Difluoro-2,4-Dinitrobenzene. Carlsberg Res. Commun. 49,
- 16. Ehrsson, H., et al. (1971). Acta. Pharm. Suecica. 8, 319.
- 17. Koshy, K.T., et al. (1975). O-(2,3,4,5,6-Pentafluorobenzyl) hydroxylamine hydrochloride as a sensitive derivatizing agent for the electron capture gas liquid chromatographic analysis of keto steroids. J. Chromatogr. Sci. 13(Feb),
- 18. Walle, T. and Ehrsson, H. (1970). Acta. Pharm. Suecica. 7, 389-406.
- 19. Benington, F., et al. (1975). Identification and separation of indolealkylamines by gas liquid chromatographic analysis of their heptafluorobutyryl derivatives. J. Chromatogr. 106, 435-439.
- 20. Kelly, R.W. and Taylor, P.L. (1976). tert-Butyldimethylsilyl ethers as derivatives for qualitative analysis of steroids and prostaglandins by gas phase methods. Anal. Chem. 48(3), 465.
- 21. Lau, H.L. (1966). J. Gas Chromatogr. 4, 136.
- 22. Tallent, W.H. and Kleinman, R. (1968). J. Lipid Res. 9, 146.
- 23. Mawhinney, T.P. and Madson, M.A. (1982). N-Methyl-N-(tert-butyldimethylsilyl) trifluoroacetamide and related N-tert-butyldimethylsilyl amides as protective silyl donors. J. Org. Chem. 47,
- 3336-3339. 24. Bazan, A.C. and Knapp, D.R. (1982). Improved derivative of 6-keto-prostaglandin  $F_{1a}$  for gas chromatographic-mass spectrometric analysis. *J. Chromatogr.* **236**, 201-207.



#### **BSTFA**

#### For excellent chromatographic separations.

$$\begin{array}{c} \text{CH}_3 \\ \text{I} \\ \text{CH}_3 - \text{Si} - \text{CH}_3 \\ \text{I} \\ \text{O} \\ \text{CH}_3 \\ \text{I} \\ \text{I} \\ \text{I} \\ \text{CF}_3 - \text{C} = \text{N} - \text{Si} - \text{CH}_3 \\ \text{CH}_3 \\ \end{array} \quad \begin{array}{c} \text{BSTFA} \\ \text{MW 257.4} \\ \text{bp 40°C/12 mm} \\ \text{d$^{20}_4$ 0.961} \\ \text{cF}_3 - \text{C} = \text{N} - \text{Si} - \text{CH}_3 \\ \text{CH}_3 \\ \end{array}$$

The greatest advantage of using Thermo Scientific BSTFA over other silvlating reagents is the increased volatility of its byproducts, mono(trimethylsilyl) trifluoroacetamide and trifluoroacetamide. This increased volatility results in the byproducts eluting with the solvent front, providing excellent chromatographic separations.

BSTFA is a powerful trimethylsilyl donor, with donor strength that is comparable to its unfluorinated analog BSA [N,O-Bis(trimethylsilyl)acetamide]. BSTFA reacts to replace labile hydrogens on a wide range of polar compounds with a -Si(CH<sub>a</sub>)<sub>a</sub> group. This physical characteristic is particularly useful in the gas chromatography of some lower boiling TMS-amino acids and TMS Krebs cycle acids.

#### **PROTOCOL**

- 1. Combine 5-10 mg sample, 0.5 ml BSTFA and 1.0 ml solvent (acetonitrile is recommended for amino acids) in a 3.0 ml Thermo Scientific Reacti-Vial\* Small Reaction Vial.
- 2. Cap vial and shake for 30 seconds.
- 3. Heat at 70°C for 15 minutes.
- 4. Analyze by gas chromatography.


NOTE: This protocol is not recommended for sugars.

Information	
Description	Pkg. Size
BSTFA [N,O-bis (trimethylsilyl) trifluoroacetamide]	25 g Hypo-Vial Sample Storage Vial
BSTFA	100 g Hypo-Vial™ Sample Storage Vial
BSTFA	10 x 1 ml ampules
	Description  BSTFA [N,O-bis (trimethylsilyl) trifluoroacetamide]  BSTFA

\* Additional hazardous handling charge.

#### BSTFA + TMCS

#### The reagent to choose for difficult-to-silylate compounds.



Thermo Scientific BSTFA + 1% TMCS is ideal for derivatizing fatty acid amides, slightly hindered hydroxyls and other difficult-to-silylate compounds. This catalyzed formulation is stronger than BSTFA alone.

#### **PROTOCOL**

- 1. Combine 5-10 mg sample, 0.5 ml BSTFA + 1% TMCS and 1.0 ml solvent (acetonitrile is recommended for amino acids) in a 3.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and shake for 30 seconds.
- 3. Heat at 7°C for 15 minutes.
- 4. Analyze by gas chromatography.

Ordering	Information	
Product #	Description	Pkg. Size
TS-38831	BSTFA + 1% TMCS [ <i>N-O</i> -bis (trimethylsilyl) trifluoroacetamide + 1% Trimethylchlorosilane]	10 x 1 ml ampules
TS-38832	BSTFA + 1% TMCS	10 g Hypo-Vial Sample Storage Vial
TS-38833	BSTFA + 1% TMCS	25 g Hypo-Vial Sample Storage Vial
<b>≭</b> TS-38834	BSTFA + 1% TMCS	100 g Hypo-Vial Sample Storage Vial
TS-38840	BSTFA + 10% TMCS [ <i>N-O</i> -bis (trimethylsilyl) trifluoroacetamide + 10% Trimethylchlorosilane]	10 x 1 ml ampules

\* Additional hazardous handling charge.

#### MSTFA and MSTFA + 1% TMCS

#### Offers maximum volatility.

#### **Highlights:**

- Trimethylsilyl donor strength comparable to BSA and
- · Reacts to replace labile hydrogens on a wide range of polar compounds with a -Si(CH3)3 group
- Used to prepare volatile and thermally stable derivatives for GC and GC/MS
- Primary advantage of Thermo Scientific MSTFA is the volatility of its byproduct, N-methyltrifluoroacetamide; MSTFA is the most volatile TMS-amide available which has an even lower retention time than MSTFA
- Often TMS derivatives of small molecules can be analyzed when derivatized with MSTFA because the byproducts and reagent itself usually elute with the solvent front
- Addition of Thermo Scientific TMCS aids derivatization of amides, secondary amines and hindered hydroxyls not derivatized by MSTFA alone

MSTFA is the most volatile TMS-amide available – its even more volatile than BSTFA or BSA. Its byproduct, N-methyltrifluoro-acetamide, has a lower retention time in GC applications than MSTFA itself. This makes it ideal for GC determinations in which the reagent or byproducts may obscure the derivative on the chromatogram. Silylation of steroids shows MSTFA to be significantly stronger in donor strength than BSTFA or BSA. MSTFA will silylate hydrochloride salts of amines directly.

#### **PROTOCOL**


- 1. Combine 5-10 mg sample, 0.5 ml MSTFA and 1.0 ml solvent (acetonitrile is recommended for amino acids) in a 3.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and shake for 30 seconds.
- 3. Heat at 70°C for 15 minutes.
- 4. Analyze by gas chromatography.

Ordering	Information	
Product #	Description	Pkg. Size
TS-48910	MSTFA ( <i>N</i> -Methyl- <i>N</i> -trimethylsilyltrifluoroacetamide)	10 x 1 ml ampules
TS-48911	MSTFA	10 g Hypo-Vial Sample Storage Vial
TS-48913	MSTFA	25 ml Hypo-Vial Sample Storage Vial
<b>≭</b> TS-48914	MSTFA	100 ml Hypo-Vial Sample Storage Vial
TS-48915	MSTFA + 1% TMCS ( <i>N</i> -Methyl- <i>N</i> -trimethylsilyltrifluoroacetamide+ 1% Trimethylchlorosilane)	10 x 1 ml ampules

★ Additional hazardous handling charge.

#### **BSA**

The perfect reagent for volatile TMS derivatives.



Under relatively mild conditions, Thermo Scientific BSA reacts quantitatively with a wide variety of compounds to form volatile, stable TMS derivatives for GC analysis. BSA is used extensively for derivatizing alcohols, amines, carboxylic acids, phenols, steroids, biogenic amines and alkaloids. It is not recommended for use with carbohydrates or very low molecular weight compounds.

BSA is used in conjunction with a solvent such as pyridine or DMF, and reactions are generally rapid. When used with DMF, BSA is the most suitable reagent for derivatizing phenols. A study of the silylating properties of BSA made by Klebe, Finkbeiner and White showed the following reactions with BSA:

- Amino acids to form both N,O-bonded TMS derivatives
- Hydroxyl compounds to form TMS ethers
- · Organic acids to form TMS esters
- Aromatic amides to form N-TMS derivatives

#### PROTOCOL 1

- 1. Combine 5-10 mg sample, 0.5 ml BSA and 1.0 ml solvent (acetonitrile is recommended for amino acids) in a 3.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and shake for 30 seconds.
- 3. Heat at 70°C for 15 minutes.
- 4. Analyze by gas chromatography.

#### PROTOCOL 2

This method was developed by E.M. Chambaz and E.C. Horning for the silylation of hydroxyl groups in sterically unhindered positions in steroids. This includes sites such as 3, 7, 16, 17(sec), 20 and 21 positions in the steroid structure. This method may be used for silylating many hydroxyl and poly-hydroxyl compounds other than steroids. It is not recommended, however, for sugars. The method is based upon the use of BSA in an uncatalyzed reaction. No trimethylchlorosilane should be used in this reaction. Hydrochlorides should be avoided because HCI also will act as a catalyst.

- 1. Combine 0.1-5.0 mg of sample and 0.2-0.4 ml BSA in a 1.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial. If material is not soluble in BSA, add 0.1-0.2 ml pyridine.
- 2. Cap vial and shake for 30 seconds.
- 3. Heat at 60°C to ease dissolution, if desired.

**NOTE:** Material is silvlated at room temperature within times varying from a few minutes to a few hours. Heating will hasten reaction.

Ordering Information			
Product #	Description	Pkg. Size	
TS-38836	BSA [N,O-bis (trimethylsilyl	)acetamide] 10 x 1 ml ampules	
TS-38838	BSA	25 g Hypo-Vial Sample Storage Vial	
× TS-38839	BSA	100 g Hypo-Vial Sample Storage Vial	

#### MTBSTFA and MTBSTFA + 1% TBDMCS

#### Offers stable TBDMS (tert-butyldimethylsilyl) derivatization.

$$\mathbf{CF_3} = \begin{matrix} \mathbf{0} & \mathbf{CH_3} & \mathbf{CH_3} \\ \mathbf{I} & \mathbf{I} & \mathbf{I} \\ \mathbf{C} - \mathbf{N} - \mathbf{Si} - \mathbf{C} - \mathbf{C} \mathbf{H_3} \\ \mathbf{I} & \mathbf{I} & \mathbf{I} \\ \mathbf{CH_3} & \mathbf{CH_3} & \mathbf{CH_3} \end{matrix} \qquad \begin{matrix} \mathbf{MTBSTFA} \\ \mathbf{MW} \ 241.3 \\ \mathbf{bp} \ 168-170^{\circ}\mathbf{C} \\ \mathbf{d}_{4}^{2} \ 1.121 \end{matrix}$$

#### **Highlights:**

- · Derivatizes hydroxyl, carboxyl, thiol and primary and secondary amines
- Typical yields are >96%
- Provides TBDMS ethers that are 104 times more stable to hydrolysis than TMS ethers
- · Reaction byproducts are neutral and volatile
- Derivatives have a high molecular concentration at M-57
- Silylating potential increased by adding 1% TBDMCS

Thermo Scientific N-Methyl-N-(tert-butyldimethylsilyl)trifluoroacetamide (MTBSTFA) provides TBDMS derivatives without the disadvantage of earlier reported TBDMS-CI formulations. Bazan and Knapp have demonstrated the usefulness of MTBSTFA by preparing an improved derivative of 6-keto-prostaglandin F1 for GC-MS.

#### **PROTOCOL**


- 1. Combine 1-10 mg of sample, 0.1 ml MTBSTFA and 0.1 ml acetonitrile in a 1.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and stand at room temperature 5-20 minutes.
- 3. Analyze by gas chromatography.

NOTE: Other solvents may be used including DMF, pyridine and THF. (DMF is not recommended for primary or secondary amines.)

Ordering Information			
Product #	Description	Pkg. Size	
TS-48920	MTBSTFA [N-Methyl-N-(tert-butyldimethylsilyl) trifluoroacetamide]	5 ml Hypo-Vial Sample Storage Vial	
TS-48927	MTBSTFA +1% TBDMCS	10 x 1 ml ampules	

#### **TMSI**

The strongest hydroxyl silylator available for carbohydrates and steroids.



Sakauchi and Horning have shown TMSI to be an all-purpose reagent for unhindered steroids to highly hindered steroids.

Thermo Scientific TMSI is unique, as it reacts guickly and smoothly with hydroxyls and carboxylic acids, but not with amines. Because TMS-derivatives are less stable than TMS-ethers or -esters, TMSI is especially useful in multiderivatization schemes for compounds containing both hydroxyl and amine groups (such as in the preparation of -0-TMS, -N-HFB derivatives of catecholamines).

TMSI is used in the derivatization of alcohols, phenols, organic acids steroids hormones glycols, nucleotides and narcotics. In addition, it is excellent for C1 through C5 fatty acids in serum and urine.

#### **PROTOCOL 1**

This method combines silylation of hydroxyl groups and acylation of amino groups. It was first used by M.G. Horning, et al. to prepare catecholamines for GC and GC/MS determinations. This method takes advantage of the fact that TMSI will silylate only hydroxyl groups. Effectively, this blocks those sites from acylation while leaving the amine sites open for acylation.

- 1. Combine and dissolve 1.0 mg sample and 1.0 ml acetonitrile in a 1.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Add 0.2 ml TMSI.
- 3. Cap vial and heat at 60°C for 3 hours.
- 4. Add 0.1 ml HFBI, TFAI or PFPI (depending on which acyl derivative is
- 5. Cap vial and heat at 60°C for 30 minutes.
- 6. Analyze by gas chromatography.

#### PROTOCOL 2

This method was developed by Sakauchi and Horning for the silylation of hydroxyl groups on highly hindered steroids. It offers fast conversion to TMS-ethers at a moderate temperature with a single reagent.

- 1. Combine 0.1-5.0 mg of steroid, 0.1-1.0 ml TMSI (0.1 ml pyridine should be added for solubilization of cortol and cortolones) in a 1.0 or 3.0 ml Reacti-Vial Small Reaction Vial.
- 2. Cap vial and heat at 100°C for 2 hours.
- 3. Analyze by gas chromatography.

#### PROTOCOL 3

- 1. Combine 400 µl TMSI and 800 µl pyridine (other solvents may be used) in a 3.0 ml Reacti-Vial Small Reaction Vial.
- 2. Add 10-15 mg sample.
- 3. Cap vial and shake until sample is dissolved. Heat to 60-70°C if needed.
- 4. Analyze by gas chromatography.

NOTE: TMSI may be used straight with carbohydrates or as a 50% solution with pyridine for wet sugars.

Ordering Product #	Information  Description	Pkg. Size
TS-88623	TMSI ( <i>N</i> -Trimethylsilylimidazole)	10 x 1 ml ampules
TS-88625	TMSI	25 g Hypo-Vial Sample Storage Vial
TS-88626	TMSI	100 g Hypo-Vial Sample Storage Vial

#### **HMDS**

The popular choice for silvlation of sugars and related substances.

Thermo Scientific HMDS greatly extends the practical range of GC, improving chromatographic results in the silvlation of sugars and related substances.

A critical study of the optimal proportions of HMDS and trimethylchlorosilane for producing maximum yield of trimethylsilyl derivatives was conducted by Sweeley, et al.

#### PROTOCOL 1

This protocol describes the method of Sweeley, et al. for the trimethylsilylation of sugars and related substances.

- 1. Combine 10 mg or less carbohydrate sample, 1.0 ml anhydrous pyridine, 0.2 ml HMDS and 0.1 ml TMCS in a 3.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and shake vigorously 30 seconds.
- 3. Let stand at room temperature 5 minutes or until derivatization is complete.
- 4. Analyze by gas chromatography.

NOTE: Solution may become cloudy when TMCS is added, due to fine precipitate of ammonium chloride. Precipitate will not interfere with gas chromatography. Carbohydrates may be warmed for 10-20 minutes at 75-85°C to hasten dissolution.

#### **PROTOCOL 2**

This method was developed primarily for silylating syrups and concentrated aqueous solutions of sugars such as starch hydrolyzates.

CAUTION: Considerable heat, ammonia gas and pressure emit during reaction. Do not premix.

- 1. Place 60-70 mg of 80% solids syrup in a 3.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Add 1.0 ml pyridine and dissolve.
- 3. Add 0.9 ml HMDS and mix.
- 4. Add 0.1 ml trifluoroacetic acid.
- 5. Shake vigorously 30 seconds.
- 6. Let stand 15 minutes.
- 7. Analyze by gas chromatography.

Ordering Information			
Product #	Description	Pkg. Size	
TS-84770	HMDS (Hexamethyldisilazane)	25 g Hypo-Vial Sample Storage Vial	
<b>≭</b> TS-84769	HMDS	100 g Hypo-Vial Sample Storage Vial	

\* Additional hazardous handling charge.

#### **TMCS**

An excellent catalyst for difficult-to-silylate compounds.

$$\begin{array}{cccc} \text{CH}_3 & \text{TMCS} \\ \text{I} & \text{MW } 108.7 \\ \text{CH}_3 - \underset{\text{CH}_3}{\text{si}} - \text{cI} & \underset{\text{d}}{\text{bp }} 57.6 ^{\circ} \text{C} \\ \text{cH}_3 & \text{d}_{\phantom{0}4}^{\phantom{0}2} \ 0.858 \end{array}$$

Thermo Scientific TMCS (trimethylchlorosilane) provides an excellent adjunct for forming trimethylsilyl ethers for GC determinations. In addition, it is used for preparing TMS derivatives of organic acids.

#### **PROTOCOL**

This protocol describes the method of Sweeley, et al. for the trimethylsilylation of sugars and related substances.

- 1. Combine 10 mg or less carbohydrate sample, 1.0 ml anhydrous pyridine, 0.2 ml HMDS and 0.1 ml TMCS in a 3.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and shake vigorously 30 seconds.
- 3. Let stand at room temperature 5 minutes or until derivatization is complete.
- 4. Analyze by gas chromatography.

NOTE: Solution may become cloudy when TMCS is added due to fine precipitate of ammonium chloride. Precipitate will not interfere with gas chromatography. Carbohydrates may be warmed for 10-20 minutes at 75-85°C to hasten dissolution.

Ordering	Information	
Product #	Description	Pkg. Size
TS-88530	TMCS (Trimethylchlorosilane)	25 g bottle

#### Methoxamine (MOX) Reagent

#### Use this reagent for preparing oximes of steroids and ketoacids prior to silylation.

Thermo Scientific MOX Reagent (M.W. 83.51) converts keto groups to methoxime derivatives. It prevents the formation of multiple derivatives (which interfere with quantitation) when enols are present during silvlation. Our MOX Reagent is a 2% solution of methoxyamine •HCI in pyridine, and it is used primarily with steroids.

The procedures below are used to prepare methoxime derivatives of steroids and ketoacids prior to silvlation. Forming methoximes is based on the work of Fates and Luukkainen, with further applications by Horning, et al. Both procedures have been used successfully by Horning, et al.

#### **PROTOCOL 1**

This simplified procedure is for stable ketones that are readily soluble in organic solvent.

- 1. Combine 2 mg sample and 0.5 ml MOX Reagent in a 10 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and heat at 60°C for 3 hours.
- 3. Add 2 ml water.
- 4. Extract with three 5 ml portions of high-purity benzene.
- 5. Combine benzene extracts and wash with 1 N HCI, followed by bicarbonate solution.
- 6. Dry over anhydrous magnesium sulfate and evaporate to 0.5 ml with
- 7. Analyze by gas or thin layer chromatography.

#### **PROTOCOL 2**

This procedure is for polar steroids, such as corticoids, that have several hydroxyl groups.

- 1. Combine 2 mg sample and 0.5 ml MOX Reagent in a 10 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and let stand overnight at room temperature.
- 3. Add 2 ml saturated NaCl solution.
- 4. Extract with three 5 ml portions of high-purity ethyl acetate.
- 5. Combine ethyl acetate extracts and wash with salt saturated 0.1 N HCI wash, followed by a 5% NaHCO3 salt saturated wash.
- 6. Dry over anhydrous magnesium sulfate and evaporate with nitrogen
- 7. Analyze by gas or thin layer chromatography.

NOTE: After completing the methoxime reaction, some researchers have silylated the reacted mixture without further treatment. The resulting mixture was centrifuged to remove solids, and aliquots of the sample were used for gas chromatography.

Ordering Information			
Product #	Description	Pkg. Size	
TS-45950	Methoxamine (MOX) Reagent (2% methoxyamine•HCI in pyridine	10 ml ) Hypo-Vial Sample Storage Vial	

#### Tri-Sil HTP (HDMS:TMCS:pyridine) Reagent

#### Our reagent-catalyst solvent mixture for one-step derivatization.

Thermo Scientific Tri-Sil HTP Reagent is composed of HMDS, TMCS and high purity pyridine. It is useful for rapid production of TMS derivatives of polar compounds for gas chromatographic determination and biochemical synthesis. The Tri-Sil HTP Reagent is ideal for GC determinations of:

- Sugars
- Alcohols
- Phenols
- Steroids
- Sterols
- · Bile acids and other organic acids
- Some amines

Our Tri-Sil HTP Reagent is based on the procedure of Sweeley, et al. and is used for the optimal conversion of organic hydroxyl and polyhydroxyl compounds into TMS ethers. The reaction proceeds as:

The versatility, speed and ease of use of our Tri-Sil HTP Reagent has made it the most widely used silylation formulation available.

#### **PROTOCOL**

- 1. Combine 5-10 mg sample and 1.0 ml Tri-Sil HTP Reagent in a 3.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Shake the reaction vigorously for 30 seconds or warm to 75-85°C to dissolve.
- 3. React at room temperature for 5 minutes.
- 4. Analyze by gas chromatography.

NOTE: A majority of hydroxyl and polyhydroxyl compounds will be completely derivatized in less than 5 minutes including sugars, phenols, organic acids, some amines and alcohols. Highly hindered compounds, such as some steroids, may require 15 minutes to 8 hours. Extremely intractable compounds may require refluxing for several hours.

Ordering Information			
Product #	Description	Pkg. Size	
TS-48999	Tri-Sil HTP Reagent HMDS:TMCS:Pyridine (2:1:10)	10 x 1 ml ampules	
<b>X</b> TS-49001	Tri-Sil HTP Reagent HMDS:TMCS:Pyridine (2:1:10)	50 ml Hypo-Vial Sample Storage Vial	

\* Additional hazardous handling charge.

#### Tri-Sil BP Reagent (BSA:pyridine)

#### A reagent-solvent solution for one-step derivatization.

Thermo Scientific Tri-Sil BP Reagent is composed of BSA (2.5 mEq/ml\*) and Pyridine.

\*1.25 mEq for amides

Tri-Sil BP Reagent reacts with:

- · Alcohols, phenols, some enols and other hydroxyl and polyhydroxyl compounds to form trimethylsilyl ethers
- Organic acids to form trimethylsilyl esters
- Aromatic amides to form N-trimethylsilyl derivatives
- Amines to form N-trimethylsilyl derivatives

In addition, Tri-Sil BP Reagent is excellent for unhindered steroids, but it is not recommended for carbohydrates.

#### **PROTOCOL**

- 1. Combine 5-10 mg sample and 1.0 ml Tri-Sil BP Reagent in a 3.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and heat at 60-70°C for 15-20 minutes to facilitate dissolution and derivatization.
- 3. Analyze by gas chromatography.

Ordering Information			
Product #	Description	Pkg. Size	
TS-49012	Tri-Sil BP Reagent (2.5 mEq/ml BSA in pyridine)	25 ml Hypo-Vial Sample Storage Vial	

#### Tri-Sil TBT Reagent (TMSI:BSA:TMCS)

#### A powerful catalyzed silylation reagent formulation.

Thermo Scientific Tri-Sil TBT is a mixture containing three parts TMSI, three parts BSA and two parts TMCS. Our Tri-Sil TBT Reagent converts all classes of hydroxyl groups to TMS ethers. Under usual conditions, the reaction is complete in a short period of time at 60-80°C. Highly hindered hydroxyls may require several hours.

#### **PROTOCOL**

This method is used to silvlate all hydroxyl groups in steroid structures, even the most sterically hindered, such as the 17 hydroxyl groups in cortol. This method also has been used by Bacon and Kokenakes to measure plasma prednisolone by GC.

- 1. Combine 0.1-5.0 mg sample and 0.2-0.4 ml Tri-Sil TBT Reagent in a 1.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and shake to dissolve.
- 3. Heat at 60-80°C for 6-24 hours to complete reaction.
- 4. Analyze by gas chromatography.

NOTE: If sample is insoluble, add 0.1-0.2 ml pyridine.

Ordering	Information	
Product #	Description	Pkg. Size
TS-49016	Tri-Sil TBT Reagent TMSI: BSA: TMCS (3:3:2)	10 x 1 ml ampule

#### Tri-Sil TP Reagent (TMSI:pyridine)

#### Great for derivatizing hydroxyl compounds.

Thermo Scientific Tri-Sil TP Reagent is a mixture of TMSI in dry pyridine (1.5 mEq/ml). It is used primarily for derivatizing hydroxyl compounds, particularly carbohydrates. Tri-Sil TP Reagent has been used successfully for the silvlation of alcohols and phenols, organic acids, hydroxylamines, amino acids, carbohydrates, flavonoids, glycols and polyglycols, nucleotides, steroids, hydroxyl acids, barbiturates, narcotics, indoles, and vitamins. Tri-Sil TP Reagent does not react with amines.

Our Tri-Sil TP Reagent can be used in the presence of water as long as there is enough reagent present to react with both the water and the sample. The reagent reads with water in a 2:1 ratio.

#### **PROTOCOL 1**

- 1. Combine 10-15 mg sample and 1.0 ml Tri-Sil TP Reagent in a 3.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and shake to dissolve. If necessary, heat at 60-70°C. Silylation is complete upon dissolution.
- 3. Analyze by gas chromatography.

#### **PROTOCOL 2**

For solutions containing ~1% or less total sugars, use a 50:50 v/v TMSI/pyridine solution.

- 1. Evaporate ~50 µl sugar solution to a glassy syrup in a 0.3 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Add 50 µl Tri-Sil TP Reagent.
- 3. Cap vial and heat at 60°C to dissolve and derivatize the sugars.
- 4. Analyze directly by gas chromatography.

Ordering Information		
Product #	Description	Pkg. Size
TS-49230	Trl-Sil TP Reagent TMSI: Pyridine (1:4)	10 x 1 ml ampules
TS-49231	Trl-Sil TP Reagent TMSI: Pyridine (1:4)	25 ml Hypo-Vial Sample Storage Vial

#### Silylation Grade Solvents

#### Manufactured to meet your exacting silylation needs.


Thermo Scientific Silylation Grade Solvents are specially manufactured and packaged to meet the exacting needs of silvlation and other sensitive derivatization reactions. Each Silvlation Grade Solvent is purified, dried and packaged under nitrogen in our convenient Hypo-Vial Sample Storage Vials. Supplied complete with elastomer septa, this unique packaging allows immediate access to your sample, without exposure to moisture and oxygen.

#### **Highlights:**

- Purified, dried and packaged under nitrogen in convenient Hypo-Vial Sample Storage Vials
- · Supplied with elastomer septa, allowing immediate access to sample without exposure to moisture and oxygen
- Use polar solvents (acetonitrile, dimethylformamide, dimethylsulfoxide, pyridine and tetrahydrofuran) to facilitate reactions; nonpolar organic solvent may be used, but they will not accelerate the rate of reaction
- Avoid water or alcohol because TMS reagents react with active hydrogen; avoid enolizable ketones
- Use dimethylformamide for steroids and other large molecules
- Use dimethylsulfoxide to prepare TMS derivatives of tertiary alcohols and some compounds with reluctant solubility in other silylation solvents
- Pyridine is an excellent solvent and reaction medium for MS reactions and is an HCI acceptor in reactions involving organochlorosilanes
- Other commonly used solvents include tetrahydrofuran and acetonitrile

Ordering Information			
Product #	Description	Pkg. Size	
<b>≭</b> TS-20062	Acetonitrile	50 ml Hypo-Vial Sample Storage Vial	
<b>≭</b> TS-20672	Dimethylformamide	50 ml Hypo-Vial Sample Storage Vial	
<b>≭</b> TS-20684	Dimethylsulfoxide	50 ml Hypo-Vial Sample Storage Vial	
<b>≭</b> TS-27530	Pyridine	50 ml Hypo-Vial Sample Storage Vial	
<b>≭</b> TS-27860	Tetrahydrofuran	50 ml Hypo-Vial Sample Storage Vial	

\* Additional hazardous handling charge. For HPLC Grade Solvents, see page 42.



Acetonitrile MW 41.05 bp 81.6°C

Dimethylformamide MW 73.09 bp 153°C

$$CH_3$$
  $S = 0$ 

**Dimethylsulfoxide** MW 78.13 bp 189°C



**Pyridine** MW 79.10 bp 115.2°C



Tetrahydrofuran MW 72.10 bp 66°C



# Servo+ and Servo

- Ready to use solutions
- High reproducibility
- High extract cleanliness
- Confidence in results

#### Acylation and Acylation Reagents

Acylation is the conversion of compounds (through the action of a carboxylic acid or a carboxylic acid derivative) that contain active hydrogens such as -OH, -SH and -NH esters; thioesters; and amides. 1 In chromatographic applications, the acylation reaction is used primarily for converting the above classes of compounds into derivatives that are better suited for chromatography<sup>2</sup> or that give a greater response to the chromatographic detection system than the parent compound.3

An important example of this application is the insertion of perfluoroacyl groups into a molecule to enhance the detectability of the substance by electron capture. The presence of a carbonyl group adjacent to the halogenated carbons enhances the electron capture detector (ECD) response.

Acyl derivatives also are useful in MS applications in which they direct the fragmentation patterns of the compounds to be studied.4

#### References

- 1. Donike, M. (1973). Acylation with bis (acylamides). N-Methyl-bis (trifluoroacetamide), two new reagents for trifluoroacetylation. J. Chromatogr. 78, 273-279.
- 2. Sullivan, J.E. and Schewe, L.R. (1977). Preparation and gas chromatography of highly volatile trifluoroacetylated carbohydrates using N-Methyl-bis (trifluoroacetamide). J. Chromatogr. Sci. 15, 196-197.
- 3. Benington, F., et al. (1975). Identification and separation of indolealkylamines by gas liquid chromatographic analysis of their heptafluorobutyryl derivatives. J. Chromatogr. 106, 435-439.
- 4. Borga, O., et al. (1971). Quantitative determination of nortriptyline and desmethylnortriptyline in human plasma by combined gas chromatography-mass spectrometry. J. Chromatogr. 4(12), 837-849.

#### Pentafluoropropanol

#### Purified for GC/MS use.



Pentafluoropropanol MW 150.05 bp 80.6°C d<sup>20</sup> 1.2880

#### **Highlights:**

- · Addition of fluorine atoms into compounds greatly enhances the sensitivity of certain detectors for those materials
- . It is advantageous to introduce fluorine atoms for ECD and GC/MS applications
- Carboxylic acids can be derivatized using a two-step reaction involving reaction with an anhydride, followed by a fluorinated alcohol

Ordering	J Information	
Product #	Description	Pkg. Size
TS-65195	Pentafluoropropanol	10 x 1 ml ampules

#### **MBTFA**

For trifluoroacylating primary and secondary amines, hydroxyl and thiol groups, and carbohydrates.

Thermo Scientific MBTFA is ideal for trifluoroacylating primary and secondary amines, hydroxyl and thiol groups, or carbohydrates under nonacidic conditions. In addition, MBTFA is used to selectively acylate amines in the presence of hydroxyl and carboxyl functions. The reaction byproduct, N-methyltrifluoroacetamide, is volatile. MBTFA also produces very volatile derivatives of carbohydrates.

#### **Highlights:**

- Trifluoroacetylates primary and secondary amines, as well as hydroxyl and thiol groups under mild nonacidic
- Principle byproduct from the derivatization reaction is N-methyltrifluoroacetamide, which is stable, volatile and does not present problems in subsequent GC
- Produces very volatile derivatives of carbohydrates and can be used to selectively acylate amines in the presence of hydroxyl and carboxyl groups that have been protected by silylation

#### PROTOCOL 1

For trifluoroacylating primary and secondary amines, and hydroxyl and thiol groups.

- 1. Combine 1-10 mg sample and 0.1-0.2 ml MBTFA. If sample is not easily solubilized, add 0.5-1.0 ml pyridine, DMF, THF or acetonitrile. (MBTFA can be pre-mixed with solvent in a 1:4 ratio. Add 1 ml pre-mixed solution to 1-10 mg compound.)
- 2. Cap vial and heat at 60-100°C for 10-30 minutes (longer for hindered compounds).
- 3. Cool and analyze by gas chromatography.

**NOTE**: MBTFA reacts with amines at room temperature to yield quantitative results in approximately 30 minutes. Hydroxyl groups are slower to react. As a general procedure, reaction mixtures should be heated for 10-30 minutes at 60-100°C. In the case of hindered compounds, further heating may be necessary.

#### **PROTOCOL 2**

For trifluoroacylating sugars.

Producing TFA derivatives of sugars using standard fluorinated anhydride and fluorinated acylimidazole procedures has yielded multiple or unstable derivatives. MBTFA produces the corresponding trifluoroacetyl derivatives of the mono-, di-, tri- and tetrasaccharides. These derivatives, when subjected to gas chromatography, produce quantitative results and yielded an unexpectedly high degree of volatility.

The high volatility of the corresponding TFA derivatives yields shorter retention times at lower temperatures than other commonly used silylation methods. The result is that polar columns with lower maximum temperature limits can be used to provide faster separations under less stringent chromatographic conditions. Mixtures of carbohydrates containing mono-through tetrasaccharides can be analyzed in a single run in as little as 15 minutes.

- 1. Combine 5-10 mg dry sugar and 0.5 ml MBTFA in a 5 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Add 0.5 ml pyridine.
- 3. Cap vial and heat at 65°C for 1 hour with occasional shaking.
- 4. Analyze 1 µl by gas chromatography.

NOTE: Reactions are complete upon dissolution.

Selective acylation of amine groups in the presence of hydroxyl and carboxyl groups is possible if these groups are first protected by silylation. The multifunctional compound first is silylated with MSTFA (N-Methyl-N[TMS] trifluoroacetamide), which silylates all of the functional groups (TMS is trimethylsilyl). The compound then is further reacted with MBTFA, exchanging the TMS-group on the amino function with a trifluoroacetyl group. Similar results are obtained with amino acids that yield N-Trifluoroacetyl-O-TMS-esters.

Ordering	Information	
Product #	Description	Pkg. Size
TS-49700	MBTFA [N-Methyl-bis(trifluoroacetamide)]	10 x 1 ml ampules
TS-49701	MBTFA	5 g bottle
TS-49703	MBTFA	25 ml
<b>≭</b> TS-49704	MBTFA	100 ml

\* Additional hazardous handling charge.

#### Perfluoroacylimidazoles (TFAI and HFBI)

Offer effective acylation of hydroxyl groups and primary and secondary amines.

R	Name	M.W.	b.p.	d <sup>20</sup> <sub>4</sub>
CF <sub>3</sub>	TFAI	164.08	38-40°C/14 mm	1.490
C <sub>3</sub> F <sub>7</sub>	HFBI	264.10	57-58°C/10 mm	1.562

In many cases, the use of N-acylimidazoles offers considerable advantages over acid chlorides and anhydrides. Advantages include:

#### **Highlights:**

- The reaction is smooth and positive, releasing no acids into the system to hydrolyze samples
- The byproduct, imidazole, is relatively inert
- Ideal for FID and ECD techniques
- Derivatives are volatile, despite size of group
- · Closely bound fluorines contribute stability

Fluorinated acylimidazoles acylate hydroxyl groups and primary and secondary amines. They react smoothly with indole alkylamines to quantitatively derivatize the indole nitrogen, as well as other functional groups present. Fluorinated imidazoles also are used for bifunctional derivatizations and in exchange reactions from the TMS derivative to the HFB derivative. In addition, a study by Ikekawa and colleagues found that O-TMS groups could be exchanged to *O*-HFB groups by adding HFBI and a small amount of HFB acid directly to the reaction mixture.

#### PROTOCOL 1

Preparing fluoro acyl derivatives for FID.

- 1. Place 0.1-2.0 mg sample in a 1.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Add 0.2 ml desired N-acylimidazole.
- 3. Cap vial and heat at 60°C for 15-30 minutes, or until reaction is complete. (Moderately and hindered steroids may require 2-6 hours heating.)
- 4. Analyze by GC using FID.

#### PROTOCOL 2

Preparing HFB derivatives of indolealkyl amines using HFBI for FID and ECD techniques.

Milligram-scale Derivatization:

- 1. Combine 2 mg sample and 0.2 ml HFBI in a 3 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and heat at 85°C for 1 hour.
- 3. Add 1 ml water and 2 ml toluene.
- 4. Cap vial and shake 2 minutes.
- 5. Analyze toluene layer by gas chromatography.

NOTE: A small amount of HFB acid remains in the toluene phase. If it interferes with analysis, wash toluene phase with 2-3 more 0.5 ml water washes.

Microgram-Scale Derivatization:

- 1. Combine 2 µg to pg quantities residue and 20 µl HFBl in a 3.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and heat at 85°C for 1 hour.
- 3. Add 2 ml pure toluene and 0.5 ml distilled water.
- 4. Cap vial and shake 2 minutes.
- 5. Remove aqueous layer.
- 6. Wash toluene layer 3 times with 0.5 ml water.
- 7. Centrifuge toluene layer 2 minutes.
- 8. Analyze toluene layer by GC using ECD.

Ordering	Information	
Product #	Description	Pkg. Size
TS-48882	TFAI (Trifluoroacetylimidazole)	10 x 1 ml ampules
<b>※ TS-44211</b>	HFBI (Heptafluorobutyrylimidazole)	5 g bottle

\* Additional dry ice and/or freight charges.

## Perfluoro Acid Anhydride (TFAA, PFAA and HFAA)

Ours are high-purity, ideal for preparing fluoracyl derivatives.

$$\begin{bmatrix} 0 & 0 \\ \parallel & \parallel \\ R - C - 0 - C - \parallel \end{bmatrix}$$

R	Name	M.W.	b.p.	$d_{\ 4}^{20}$
CF <sub>3</sub>	TFAA	210.0	39.5-40.5°C	1.490
C <sub>2</sub> F <sub>5</sub>	PFAA	310.0	74°C	1.571
C <sub>3</sub> F <sub>7</sub>	HFAA	410.0	106-107°C	1.665

Fluorinated anhydrides are used to prepare fluoracyl derivatives for GC/MS, they produce stable volatile derivatives for FID and ECD techniques.

#### **PROTOCOL 1**

Preparing fluoracyl derivatives of amines and alcoholic compounds on a submicrogram scale for ECD.

- 1. Combine < 50 ng sample dissolved in 500  $\mu l$  benzene and 100  $\mu l$  0.05 M TEA in benzene in a 5.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Add 10 µl Acid Anhydride.
- 3. Cap vial, heat at 50°C for 15 minutes, then cool.
- 4. Add 1 ml water, cap vial and shake 1 minute.
- 5. Add 1 ml 5% aqueous ammonia, cap vial and shake 5 minutes.
- 6. Centrifuge.
- 7. Inject 1-10 µl benzene phase for ECD.

NOTE: Use 250 µg for FID. Excess TEA is required for quantitative conversion of amines. TEA does cause disturbances in the chromatogram at high EC sensitivity. The benzene is used as sample solvent and TEA solvent should be dry, as water will compete for the anhydride during reactions. The amount of anhydride used in the procedure (10  $\mu$ I) is 25% more than necessary for a complete reaction - even if the 0.5 ml benzene used is water-saturated.

#### **PROTOCOL 2**

Preparing fluoracyl derivatives of phenols for FID and ECD.

For Flame Ionization Detection:

- 1. Combine 1 mg sample dissolved in 0.5 ml benzene and 200  $\mu$ l 0.1 M TEA in
- 2. Add 25 µl Acid Anhydride.
- 3. Cap vial and let react at room temperature for 15 minutes.
- 4. Add 0.5 ml 1 M phosphate buffer, pH 6.0, and shake for 30 seconds.
- 6. Separate organic phase. Analyze by GC.

For Electron Capture Detection:

- 1. Combine 0.5 ml benzene containing the sample and 10 µl TEA in benzene in a 5 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Add 10 µl HFAA.
- 3. Cap vial and let react at room temperature for 10 minutes.
- 4. Add 0.5 ml 1 M phosphate buffer, pH 6.0, and shake 30 seconds.
- 5. Centrifuge; analyze 2 µl benzene phase by GC.

NOTE: Excess anhydride and acid are removed by the aqueous extraction. No ECD disturbances from the water or other constituents. HFB-esters of phenols are stable in the presence of water (with aqueous phase at  $pH \le 6.0$ ). Alkaline extraction with reagents, such as aqueous ammonia, decomposes the HFB ester. A pH  $\leq$  6.0 maintains the TEA catalyst in the protonized form. TEA in the unprotonized form will catalyze decomposition of the esters.

Ordering	Information	
Product #	Description	Pkg. Size
<b>≭</b> TS-67363	TFAA (Trifluoroacetic Acid Anhydride)	100 g bottle
TS-65193	PFAA (Pentafluoropropionic Acid Anhydride)	10 x 1 ml ampules
<b>≭</b> TS-65192	PFAA	25 g bottle
<b>X</b> TS-65191	PFAA	100 g bottle
TS-63164	HFAA (Heptafluorobutyric Acid Anhydride)	10 x 1 ml ampules
<b>≭</b> TS-63163	HFAA	25 g bottle
<b>≭</b> TS-63162	HFAA	100 g bottle

\* Additional hazardous handling charge.

#### Alkylation and Alkylation Reagents

When used in derivatization for gas chromatography, alkylation represents the substitution of an active hydrogen by an aliphatic or aliphatic-aromatic1 (benzyl) group. This technique is used to modify those compounds containing acidic hydrogens, such as carboxylic acids and phenols. The principal chromatographic use of this reaction is the conversion of organic acids into esters, which produce better chromatograms than the free acids.

In addition, alkylation reactions can be used to prepare ethers, thioethers and thioesters; N-alkylamines; and amides.<sup>2</sup> As the acidity of the active hydrogen decreases. the strength of the alkylating reagent must be increased. As the reagents and conditions become harsher, the selectivity and applicability of the methods become more limited.

#### References

- 1. Kawahara, F.K. (1968). Microdetermination of derivatives of phenols and mercaptans by means of electron capture gas chromatography. Anal. Chem. 40(6), 1009.
- 2. Kananen, G., *et al.* (1972). Barbiturate analysis a current assessment. J. Chrom. Sci. 10, 283-287.

# BF<sub>3</sub>-Methanol

Provides convenient, fast and quantitative esterification of fatty acids.

BF<sub>a</sub>-Methanol methylation is one of the most convenient ways to prepare methyl esters of fatty acids. Supplied in an easy-to-use septum-sealed Hypo-Vial Sample Storage Vial, Thermo Scientific BF<sub>3</sub>-Methanol offers convenient syringe removal of your sample – without exposing the contents.

#### **PROTOCOL 1**

For preparing fatty acid methyl esters.

- 1. Combine 100 mg fatty acid and 3 ml BF<sub>3</sub>-Methanol in a 5.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and heat at 60°C for 5-9 minutes.
- 3. Cool and transfer to separatory funnel with 30 ml hexane.
- 4. Wash 2 times with saturated NaCl solution.
- 5. Discard aqueous layers.
- 6. Dry over sodium sulfate.
- 7. Evaporate solvent under nitrogen.
- 8. Analyze by gas chromatography.

#### PROTOCOL 2

For preparing C8-C17 fatty acids.

- 1. Combine 500 mg fatty acid and 5 ml BF<sub>2</sub>-Methanol in a 25 ml flask.
- 2. Heat on a steam bath for 5 minutes.
- 3. Add saturated NaCl solution until total volume is ~20 ml.
- 4. Cap flask and invert several times.
- 5. Allow organic layer to collect at the top, then separate.
- 6. Dry organic layer over Na<sub>2</sub>SO<sub>4</sub>.
- 7. Evaporate solvent under nitrogen.
- 8. Analyze by gas chromatography.

Ordering	Information	
Product #	Description	Pkg. Size
<b>≭</b> TS-49370	BF <sub>3</sub> -Methanol	100 ml Hypo-Vial Sample Storage Vial

\* Additional hazardous handling charge.

#### Pentafluorobenzyl Bromide (PFBBr)

For the electron capture GC analysis of carboxyl acids, phenols and sulfonamides.

Pentafluorobenzylation by an "Extraction Alkylation" technique has been described for the electron capture GC analysis of carboxyl acids, phenols and sulfonamides. This process uses tetrabutylammonium as a counter ion and methylene chloride as a solvent. Reaction times are fast (~20 minutes). Derivatives are highly EC-sensitive, making them useful in low-level determinations of fatty acids.

Kawahara performed extensive work with this reagent, using a potassium carbonate catalyst for the electron capture analysis of mercaptans, phenols and organic acids in surface water.

PFBBr has been applied in analyzing trace organics in asphalts, as a "fingerprinting" technique for identifying asphalt pollutants found in surface waters.

#### **PROTOCOL**

For pentafluorobenzylation of acids, phenols and sulfonamides.

- 1. Place 1 ml methylene chloride containing 0.2 mg sample in a 3 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Add 1 ml aqueous 0.1 M TBA hydrogen sulfate and 0.2 M sodium hydroxide solution.
- 3. Add 20 µl PFBBr.
- 4. Cap vial and shake for 20-30 minutes.
- 5. Inject a portion of the methylene chloride phase into the chromatograph for FID analysis
- 6. Evaporate methylene chloride to dryness with nitrogen and redissolve with benzene for ECD analysis.

Ordering	Information	
Product #	Description	Pkg. Size
TS-58220	PFBBr (Pentafluorobenzyl Bromide)	5 g

#### Methylate Reagent (DMFDMA)

For easy, effective preparation of methyl esters from fatty acids and amino acids.

For preparing methyl esters for gas chromatography, Thermo Scientific Methylate Reagent offers significant advantages including:

- Speed the reaction is complete upon dissolution
- Quantitation quantitative yields are obtained when reagent and sample are injected - without prior mixing
- Your choice of formulation our Methylate Reagent is a convenient, ready-to-use reagent that contains 2 mEq/ml in pyridine.

Our Methylate Reagent is stable at room temperature and is packed in convenient, ready-to-use Hypo-Vial Sample Storage Vials. No water washing, extraction or concentration of the derivatives are required. Plus, no water is formed in the reaction.

Reactions with the Methylate Reagent usually are complete upon dissolution, except for long chain solid acids. In these applications, it is necessary to use Methylate Reagent with additional solvent and mild heating. Suitable solvents include pyridine, benzene, methanol, chloroform, methylene chloride, THF and DMF.

Thenot, et al. have demonstrated analytical applications that use of the Methylate Reagent for analyzing fatty acids1 and amino acids.2

#### **PROTOCOL 1**

Methods of alkylation using DMF-Dialkyl Acetal Reagents.

- 1. Combine 50 mg fatty acid and 1 ml Methylate Reagent in a 3 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and heat at 60°C for 10-15 minutes or until dissolution is complete.
- 3. Analyze by gas chromatography.

#### **PROTOCOL 2**

For preparing N-dimethylaminomethylene (DMAM) methyl esters of

- 1. Combine amino acid with 1:1 ratio of Methylate Reagent to acetonitrile in a Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Cap vial and heat at 100°C for 20 minutes or until dissolution is complete.
- 3. Analyze by gas chromatography.

NOTE: Aspartic acid requires a longer time for complete dissolution. Hydroxyl groups on hydroxyl-substituted amino acids do not react under the above conditions.

Product #	Description	Pkg. Size
TS-49350	Methylate Reagent (2 mEq/ml in pyridine) (N,N-Dimethylformamide dimethyl acetal)	25 ml Hypo-Vial Sample Storage Vial

#### MethElute Reagent (TMPAH)

A powerful reagent for accurate, sensitive on-column methylation.

$$\begin{array}{c} \text{CH}_3 \\ \text{[CH}_3-\text{N-CH}_3]^{\oplus} \end{array} \begin{array}{c} \ominus \\ \text{[OH]} \end{array} \begin{array}{c} \bigoplus \\ \text{heat} \end{array} \begin{array}{c} \text{[CH}_3]^{\oplus} \\ \text{Ion} \end{array} \begin{array}{c} \bigoplus \\ \text{Methylcarbonium} \\ \text{Ion} \end{array} \begin{array}{c} \bigoplus \\ \text{Methylcarbonium} \\ \text{Methylated} \\ \text{Phenobarbital} \end{array}$$

Figure 1. Thermo Scientific MethElute Reagent reaction with phenobarbital.

Thermo Scientific MethElute Reagent [0.2 M trimethylanilinium hydroxide (TMPAH) in methanol solution] is a powerful methylating reagent for quantitative methylation and detection of barbiturates, sedatives, xanthine bases, phenolic alkaloids and Dilantin by gas chromatography.

Our MethElute Reagent gives a single quantitative peak response for each substance. When the reagent is heated with drug-containing extracts, serum or urine, those drugs containing reactive amino, hydroxyl and carboxyl functions will be methylated at the reactive sites.


#### **Performance Characteristics**

Accuracy. Comparable or better than the UV/TLC method. When phenobarbital and Dilantin are present, the UV/TLC method cannot be used, as Dilantin interferes with the phenobarbital determination. The GC procedure yields excellent results for this combination of drugs.

Precision. The coefficient of variation is 5% or less.

Sensitivity. Detects barbiturates down to 0.2 mg/dl.

Ordering Information		
Product #	Description	Pkg. Size
TS-49300	MethElute Reagent (TMPAH)	10 ml Hypo-Vial Sample Storage Vial
TS-49301	MethElute Reagent (TMPAH)	12 x 1 ml Hypo-Vial Sample Storage Vials



# Redefine your expectations

- High temperature stability
- Low bleed
- Superior inertness
- Improved reproducibility

# TraceGOLD GC Columns www.thermoscientific.com/tracegold

#### Siliconizing Fluids

#### An easy way to bond polymer films to surfaces.

Thermo Scientific Multifunctional Siliconizing Fluids are specially designed to chemically bind microscopically thin, water-repellent films to glass, quartz, silica and ceramics. The coated surfaces are neutral, hydrophobic and non-oily. In addition, they offer increased resistivity and are not affected by solvents that are not readily hydrolyzed.

#### Use our Siliconizing Fluids to treat pipettes, beakers, certain plastics, ceramics, fiber optics and more:

- For clean drainage and elimination of meniscus
- To reduce adsorption of polar compounds, proteins and trace metals onto glass surfaces and reduce leaching of trace metals into solution
- To prevent current tracking and minimize electrical leakage on glass surfaces and ceramics
- To protect delicate samples against the possible reactive effects of -OH sites present on all types of glassware

#### More reasons to use our Siliconizing Fluids:

- Easy to apply
- Economical a little bit goes a long way
- Surfaces can easily be recoated
- Improve surface water-repellency
- Increase surface resistivity

#### Siliconizing Fluid – Water Soluble

Our water-soluble fluid for siliconizing glass surfaces.

## $\mathrm{CH_3(CH_2)_{16}CH_2Si(OR)}$ Siliconizing Fluid – Water Soluble

Thermo Scientific Water Soluble Siliconizing Fluid is an easy-to-use silane monomer solution that is supplied as a 20% solid solution in a mixture of diacetone alcohol and tertiary butyl alcohol. The primary silane component is an octadecyltrialkoxysilane that, when mixed with water, is hydrolyzed to a silanol. This silanol condenses with available hydroxyl groups and other silane monomers to form a film on the glass surface.

Our Water Soluble Siliconizing Fluid is especially useful in the biochemical field because of its aqueous phase application to glass and its greater resistance to base hydrolysis.

#### **Instructions for Use**

Our Water Soluble Siliconizing Fluid is applied to a clean glass surface as a dilute aqueous solution. Prepare a 0.1-1.0% solution by weight or volume. [Note: Siliconzing Fluid - Water Soluble contains 20% solids; therefore, one part our Water Soluble Siliconizing Fluid plus 99 parts water (w/w) yields a 0.2% solution - not 1%.] Add the Siliconizing Fluid to water with constant stirring. A clear-to-slightly hazy solution will be obtained. Solutions are not indefinitely stable and will turn cloudy and precipitate after several days and, therefore, solution should be prepared just before use. Solution stability can be optimized, however, by adjusting the aqueous solution pH to 4.5-5.0.

The article to be coated is dipped into the solution, or the surface is flooded with the solution. A thin film will immediately lubricate the glass surface, making it water-repellent. The surface then is air-dried for 24 hours or heated at 100°C for several minutes. Exact drying conditions should be determined before use in commercial process applications.

Ordering	Information	
Product #	Description	Pkg. Size
<b>≭</b> TS-42799	Siliconizing Fluid – Water Soluble	120 ml

X Additional hazardous handling charge.

#### Siliconizing Fluid - Hydrocarbon Soluble

#### A hydrocarbon soluble fluid for siliconizing glass surfaces.

Siliconizing Fluid -**Hydrocarbon Soluble** 

Thermo Scientific Hydrocarbon Soluble Siliconizing Fluid is a short chain, clear polymeric silicone fluid consisting primarily of dichlorooctamethyltetrasiloxane. When applied to glass, guartz or similar products, the unhydrolyzed chlorines present on the chain react with surface silanols to form a neutral, hydrophobic and tightly bonded film over the entire surface.

Our Hydrocarbon Soluble Siliconizing Fluid is ideal for use on metals, certain plastics, ceramics and fiber optics.

Our Hydrocarbon Soluble Siliconizing Fluid is acidic and care should be taken to avoid corrosion of metal that comes into contact with the liquid. The fluid is acidic only during application. After application the surface is neutral.

Caution: Material is flammable before film is formed and HCI fumes are generated in the application.

#### **Instructions for Use**

Wipe-on treatment: Wearing rubber gloves, wet a cloth with undiluted fluid and rub it on the clean surface until an oily film is formed. Rub with a dry cloth until the surface is clear.

Solution treatment: Dilute Hydrocarbon Soluble Siliconizing Fluid with 1-10% clean dry solvent such as acetone, toluene, carbon tetrachloride, methylene chloride or hexane. Do not use esters or alcohols. Articles can be dipped and air-dried. No heating is required. For a slightly more durable coating, heat articles at 100°C for 30 minutes.

Our Hydrocarbon Soluble Siliconizing Fluid (as supplied) is stable for at least one year. Discard prepared solutions after use.

Jevons, et al. reported treating all glassware with a 10% v/v solution of Siliconizing Fluid - Hydrocarbon Soluble in carbon tetrachloride.1

Ordering Information		
Product #	Description	Pkg. Size
<b>≭</b> TS-42800	Siliconizing Fluid-Hydrocarbon Soluble	120 ml
<b>≭</b> TS-42801	Siliconizing Fluid-Hydrocarbon Soluble	480 ml
<b>≭</b> TS-42855	Siliconizing Fluid-Hydrocarbon Soluble	5 x 10 ml ampules

\* Additional hazardous handling charge.

#### Monofunctional Silane - HMDS

#### Ideal for active surface site deactivation.

$$\begin{array}{c|cccc} \text{CH}_3 & \text{CH}_3 & \text{HMDS} \\ \text{I} & \text{I} & \text{MW 161.4} \\ \text{CH}_3 & \text{Si} - \text{NH} - \text{Si} - \text{CH}_3 & \text{MW 161.4} \\ \text{I} & \text{I} & \text{bp 125°C} \\ \text{CH}_3 & \text{CH}_3 & \text{n}_3^{20} \text{1.4071} \end{array}$$

Thermo Scientific HMDS is a popular monofunctional silane that many researchers have found useful for deactivating and coating chromatographic supports. Because of its monofunctional nature, this silane can react with only one site on the surface. Polymerization is not possible, eliminating the chances for unbound polymers to float free and elute from the column - avoiding exposure of unreacted silanols beneath the layer. In addition, surface moisture is eliminated because monofunctional reagents dehydrate the surface.

There are several methods for deactivating surfaces with HMDS:

- 1. Slurrying or dipping the items to be deactivated in a 5-10% solution of the reagent in an unreactive solvent.
- 2. Vapor phase deactivating by pulling straight vapor into an evacuated container that holds the item to be deactivated.
- 3. Placing the item and a few milliliters of the reagent in a beaker, then placing a watch glass on top (as in the case of glass wool silanization).

#### References

- 1. Nawrocki, J. (1985). Modification of silica with mixture,
- at hexamethylcyclotrisilazane. Chromatographia 20(5). 2. Owens, N.F., et al. (1987). Inhibition of cell adhesion by a synthetic polymer
- absorbed to glass shown under defined hydrodynamic stress. J. Cell Sci. 87, 667-675.

Ordering Information		
Product #	Description	Pkg. Size
<b>≭</b> TS-84769	HMDS (Hexamethyldisilazane)	100 g
TS-84770	HMDS (Hexamethyldisilazane)	25 g

\* Additional hazardous handling charge

# **Chromatography Resources**

Our 2012-2013 Chromatography Columns and Consumables catalog helps you find all your chromatography needs in a single resource.

To order your copy visit www.thermoscientific.com/catalog

Our bi-monthly Separated by Experience newsletter keeps you up-to-date on the latest technical and product information of interest to chromatographers.

Subscribe today! www.thermoscientific.com/chromatography



#### Handbook of Analytical Derivatization Reaction

A self-contained methodology reference manual and efficient entry point to the original literature resource book.



The Handbook of Analytical Derivatization Reactions by Daniel R. Knapp is a general collection of analytical derivatization methods for chromatography and mass spectroscopy involving the formation of covalent derivatives before analysis.

Methods contained in this volume are organized according to the type of sample being derivatized. Methods include structural formulas, experimental directions and useful comments. A thorough system of indexing takes you quickly to the "lab ready" methods of interest.

#### Product # Description

TS-15012 Handbook of Analytical Derivatization Reactions Knapp, D.R. Ed (1979) Published by John Wiley and Sons, Inc. Hardcover, 741 pages



#### Introduction to HPLC Ion Pair Reagents

High-purity reagents with the selectivity needed for good separation.

In the past, reverse-phase HPLC analysis of highly charged acidic and basic compounds was frustrating and resulted in poor resolution. Important biomolecules such as amino acids, peptides, organic acids, polyamines and catecholamines had to be separated by ion exchange or by suppression techniques.

Thermo Scientific Ion Pair Reagents enable you to quickly and efficiently analyze charged compounds using reverse-phase techniques. Our ion pair reagents are simply dissolved in the HPLC solvent system, resulting in the formation of stable chromatographic complexes that can be separated using reverse-phase columns. By using the correct ion pair reagents, you achieve:

- Increased or decreased retention, permitting controlled selectivity
- Resolution of complex ionic mixtures without using ion exchange columns
- · Improved peak symmetry

#### Reverse-phase ion pair chromatography theories

Two principal theories have been proposed to explain reverse-phase ion pair chromatography. In the first theory, small polar ion pair reagents react with the ionized solute, forming neutral ion pairs. The second theorizes that an active ion exchange surface is produced in which long chain, nonpolar anions and cations are absorbed by the hydrophobic stationary phase.

To optimize chromatographic separations in ion pair elution systems, high-purity reagents of exceptional optical transparency are needed. Ion Pair Reagents are specially purified for ion pair chromatography and provide the selectivity needed for good separations.

- 1. Bennett, H.P.S., et al. (1981). Biochemistry. 20, 4530.
- 2. Starratt, A.N. and Stevens, M.E. (1980). J. Chromatogr. 194, 421.
- 3. Burgess, A.W., et al. (1982). Proc. Natl. Acad. Sci. U.S.A. 79, 5753.
- 4. Hearn, M.T.W. and Grego, B. (1983). J. Chromatogr. 255, 125.
- 5. Shoneshofer, M. and Fenner, A. (1981). J. Chromatogr. 224, 472.
- 6. Fischli, W., et al. (1982). Proc. Natl. Acad. Sci. U.S.A. 79, 5435.
- 7. Hancock, W.S., et al. (1979). J. Chromatogr. **168**, 377. 8. Hearn, M.T.W., et al. (1978). J. Chromatogr. **157**, 337.
- 9. Hearn, M.T.W. and Grego, B. (1983). *J. Chromatogr.* **255,** 125.
- 10. Hearn, M.T.W. and Grego, B. (1983). J. Chromatogr. 266, 75.
- 11. Rivier, J. (1978). J. Liq. Chrom. 1, 343.

#### Heptafluorobutyric Acid

An ion pair reagent for the reverse-phase HPLC separation of proteins and peptides.

#### **Highlights:**

- · Clear, colorless liquid
- Typical purity is 99.7% by GC; less than 0.1% water
- Sequencing reagent for classical and automated Edman degradation of peptides and proteins
- Density: 1.645; B.P.: 120°C
- · Packaged under nitrogen in amber glass ampules or bottles

#### References

- 1. Hearn, M.T.W. and Hancock, W.S. (1979). Trends Biochem. Sci. 4, N58-N62.
- 2. Bennett, H.P.J., et al. (1980). J. Liquid Chromatogr. 3, 1353-1366.
- 3. Bennett, H.P., et al. (1981). Biochemistry 20, 4530-4538.

Ordering Information		
Product #	Description	Pkg. Size
<b>≭</b> TS-25003	Heptafluorobutyric Acid, Sequencing Grade	100 ml
TS-53104	Heptafluorobutyric Acid, HPLC Grade	10 x 1 ml

X Additional hazardous handling charge.

#### Triethylamine (TEA)

#### Ideal for HPLC separation and analysis of peptides!

Triethylamine is an ion-pairing reagent that alters selectivity in reverse-phase HPLC separations. By pairing with peptides, it effectively sharpens peaks, resulting in improved peak resolution. Triethylamine comes in two grades. Our ionate grade is designed for use as an ion-pair reagent in HPLC separations and has a low UV absorbance to provide you with the most sensitive detection across all wavelengths. Seguanal grade is designed to meet the special requirements for peptide sequencing and analysis.

#### **Highlights:**

- > 99.5% triethylamine purity, allowing sensitive peptide detection at low UV wavelengths in reverse-phase HPLC peptide separation systems
- Packaged in amber glass bottles with protective Teflon® TFE-lined fluorocarbon caps for reagent integrity




Figure 1A-1D. Effect of TEA concentration on a mixture of basic antihistamines and 2,4-dimethylaniline\* 150 mm x 4.6mm C8. Conditions: a) 40% methanol, 0.060 M HSA sodium salt, 0.045 M citric acid; b) 0.150 M citric acid, 0.060 M TEA, pH 7.5 with NaOH; c) 0.150 M citric acid, pH 7.5 with NaOH; isocratic with TEA concentrations modified by varying b/c ratio, 3 min./ml, 50°C, 254 nm.

Ordering Information		
Product #	Description	Pkg. Size
TS-53101	Triethylamine (TEA), HPLC Grade	25 g
<b>≭</b> TS-25108	Triethylamine (TEA), Sequencing Grade	100 g

\* Additional hazardous handling charge.

#### Formic Acid Ampules

#### Ideal reagent for LC-MS applications.

High Purity Thermo Scientific Formic Acid is sealed in amber glass ampules under a dry nitrogen atmosphere. A pre-measured aliquot of acid greatly simplifies the preparation of liter quantities of mobile phases at the standard 0.1% formic acid concentration. The quality of the formic acid coupled with ampule packaging provides reliability and convenience that adds value to both the chromatographic and MS results.

#### Formic Acid, Reverse-Phase HPLC and Mass Spectrometry

Formic acid is a component commonly found in reverse-phase mobile phases to provide protons for LC/MS analysis. The presence of a low concentration of formic acid in the mobile phase is also known to improve the peak shapes of the resulting separation. Unlike trifluoroacetic acid (TFA), formic acid is not an ion-pairing agent and it does not suppress MS ionization of polypeptides when used as a mobile-phase component.

#### **Highlights:**

#### 99% pure formic acid

- Consistent LC baselines
- No potential interference introduced in LC or MS applications
- · No signal suppression in the mass spectrometer

**Formic Acid** MW 46.03

#### High-performance ampule packaging

· Amber glass, pre-scored, nitrogen-flushed ampules protect formic acid from light and moisture

#### **Convenient format**

 Ampule packaging simplifies the preparation of gradient and isocratic mobile phases containing 0.1% (v/v) formic acid in water or acetonitrile; the contents of a single vial in a final volume of 1 L of solvent yields a mobile phase of the most common formic acid concentration.

Ordering	Information	
Product #	Description	Pkg. Size
TS-28905	Formic Acid 99+%	10 x 1 ml ampules

#### Trifluoroacetic Acid (TFA)

#### 1 ml ampules allow you to make a fresh 0.1% TFA solution in seconds!

Thermo Scientific Trifluoroacetic Acid (TFA) is the most commonly used ion pairing agent in reverse-phase peptide separations because TFA:

- · Sharpens peaks and improves resolution
- Is volatile and easily removed
- · Has low absorption within detection wavelengths
- Has a proven history

#### **Highlights:**

#### **Purity**

• > 99.5% TFA purity and exceptional clarity, allowing sensitive, nondestructive peptide detection at low UV wavelengths in reverse-phase HPLC protein and peptide separation systems

#### **High-performance packaging**

• Our TFA is packaged under nitrogen in amber glass ampules or bottles with protective Teflon - TFE-lined fluorocarbon caps to ensure TFA integrity

#### **Economical convenience**

 Choose the TFA format that works best for your application. In just a few seconds, 1 ml ampules can be used to prepare 1 liter of fresh 0.1% v/v trifluoroacetic acid solution for the mobile phase in reverse-phase chromatography

#### **Applications:**

- Ion pair reagent for reverse-phase HPLC1-3
- Protein/peptide sequencing<sup>4-7</sup>
- Protein/peptide solubilizing agent<sup>4-7</sup>
- Solid-phase peptide synthesis<sup>8</sup>
- · Amino acid analysis

#### **Making 0.1% Solutions of Trifluoroacetic Acid**

For complex peptide separations, the key to success can be to vary selectivity. Varying mobile phase composition on the same column can change selectivity enough to resolve peptides that would otherwise overlap. Trifluoroacetic acid is the most frequently used modifier for peptide separations in reverse-phase HPLC. The TFA concentration usually specified is 0.1%. For reproducible separations from run-to-run or from lab-to-lab, it is essential to make TFA concentrations the same.

Trifluoroacetic acid concentration can and should be specified as either "w/v" (weight/volume), or as "v/v" (volume/volume). The w/v specification designates that the TFA is to be weighed and added to a volume of mobile phase (e.g. 0.1% TFA w/v requires one gram of TFA per liter). The v/v specification designates that the TFA is to be measured by volume (e.g. 0.1% TFA v/v requires one ml of TFA per liter).

Because the density of trifluoroacetic acid is 1.53 g/ml the difference between 0.1% TFA (w/v) and 0.1% TFA (v/v) is more than 50%. For the sake of reproducibility, it is essential for authors of a method to specify, and for users of a method to know, whether the TFA concentration is given as "w/v" or "v/v".


## **Thermo Scientific Hypercarb Columns**

100% porous graphitic carbon for extended separation capabilities


- Exceptional retention of very polar analytes
- · Separates structurally related substances
- pH stable from 0 to 14
- Ideal for high temperature applications

#### www.thermoscientific.com/hypercarb






A planar compound can align itself closely with the Hypercarb® surface resulting in more interaction and retention




A non-planar molecule cannot align itself as closely with the Hypercarb surface resulting in less interaction and retention

Electrostatic charge of polar analyte approaching the graphite surface



Electrons in the polarizable surface of the graphite attracted towards positive charge



Electrons in the polarizable surface of the graphite repelled away from a negative charge

#### References

- 1. Chicz, R.M. and Regnier, F.E. (1990) Methods Enzymol. 182, 392-421.
- 2. Zarember, K.A., et al. (2002) Infect. Immun. 70, 569-576.
- 3. Lassy, R.A. and Miller, C.G. (2000) J. Bacteriol. 182, 2536-2543.
- 4. Smith, B.J. (1997) Protein Sequencing Protocols. Humana Press.
- 5. Allen, G. (1989) Sequencing of Proteins and Peptides, Second Revised Addition.
- 6. Backstrom, J.R., et al. (1996) J. Neurosci. 16, 7910-7919.
- 7. Hermann, P.M., et al (2000) J. Neurosci. 20, 6355-6364.
- 8. Stuart, J.M. and Young, J.D. (1984) Solid Phase Peptide Synthesis, Second Edition. Pierce Chemical Company.

Ordering	Information	
Product #	Description	Pkg. Size
<b>≭</b> TS-28901	Trifluoroacetic Acid, Sequencing Grade	500 ml
<b>≭</b> TS-28902	Trifluoroacetic Acid, Sequencing Grade	10 x 1 g
<b>≭</b> TS-28903	Trifluoroacetic Acid, Sequencing Grade	100 g
<b>≭</b> TS-28904	Trifluoroacetic Acid, Sequencing Grade	10 x 1 ml ampules

<sup>\*</sup> Additional hazardous handling charge.

#### Thermo Scientific Derivatization Reagents for HPLC

#### Designed to provide selectivity and improve sensitivity.

The lack of a universal HPLC detector that provides high sensitivity (as well as some degree of selectivity) established the need for suitable derivatization procedures. Derivatization is the chemical modification of an existing compound, producing a new compound that has properties more suitable for a specific analytical procedure. It is an analytical tool that can be used to provide both selectivity and improved sensitivity.

There are several requirements for a derivatization protocol:

- 1. At least one acidic, polar functional group must be available for reaction on the parent compound.
- 2. A single derivative should be formed per parent compound.
- 3. The reaction should be reproducible under the given time and reaction conditions.
- 4. The reaction should proceed quickly and easily under mild conditions.
- 5. The reaction byproducts (if any) should not interfere with the chromatography, or with detection of the sample.

Pre- and post-chromatographic techniques are both used in HPLC derivatization. In addition, off-line and on-line reactions have been employed with both techniques. Pre-chromatographic (or pre-column techniques) offer more than greater selectivity and sensitivity in detection. Pre-column techniques can be used to enhance stability, improve resolution, improve peak symmetry and increase or decrease retention of solutes. FDAA (Marfey's Reagent) allows separation and quantification of optical isomers of amino acids (Figure 2). Post-chromatographic (or post-column) techniques are used primarily to provide selectivity and improve sensitivity.

We offer a variety of HPLC detection reagents for pre- and post-chromatographic techniques. All compounds and formulations are purified for chromatography, minimizing artifact formation.




Figure 1. Separation of D- and L-DOPA on 100 mm x 4.6 mm C18. Conditions: A) 0.05 M triethylamine phosphate, pH 3.0; B) acetonitrile. Linear gradient: 10 to 40% B in 45 minutes, 2.0 ml/minute, 25°C, 340 nm.

Thermo Scientific	Detection Reagents for HPLC			
Functional Group	Description	Detection*	Page	Comments
Carboxylic Acid 0 II R – C – OH	<i>p</i> -Bromophenacylate  O  Br C - CH <sub>2</sub> -Br	UV	31	Formulation of 1.0 mmol/ml $p$ -bromophenacyl bromide and 0.005 mmol/ml crown ether in acetonitrile; pre-column; nanomole detection levels: $\lambda_{max} = 260 \text{ nm}^{1.7}$
Primary Amine R-N-H H	Dabsyl Chloride  N  N  S  CI	Vis	41	4-N, N-dimethylaminoazobenzene-4'-sulfonyl chloride (dabsyl chloride); pre-column; nanomole detection levels: $\lambda_{\text{max}} = 436 \text{ nm}^{\text{8-14}}$
	FDAA, Marfey's Reagent  ON ON HOND NH2	UV	31, 41	1-fluoro-2,4-dinitrophenyl-5-L-alanine amide (FDAA); pre-column; nanomole detection levels: $\lambda_{\rm max}=340$ nm. For chiral separations of amino acids. <sup>15, 28-23</sup>
	Ninhydrin OHOHOH	Vis	37	Post-column; nanomole detection levels: $\lambda_{\text{max}}$ = 570 nm <sup>22</sup>
	PITC No Constant N	UV	40	Phenylisothiocyanate (PITC); pre-column; picomole detection levels: $\lambda_{\text{max}} = 254 \text{ nm}^{23\text{-}24}$
	TNBSA  0   0   0   0   0   0   0   0   0   0	EC, UV	35	2,4,6-Trinitrobenzene-sulfonic acid (TNBSA); pre- or post-column; nanomole detection levels with EC and UV, GC - 0.85V; $\lambda_{max}=250~\text{nm}^{25-26}$
Secondary Amine R—NH—R'	Ninhydrin (see structure above)	Vis	37	Post-column; nanomole detection levels: $\lambda_{\rm max}$ = 440 nm <sup>22</sup>
	PITC (see structure above)	UV	40	Phenylisothiocyanate (PITC); pre-column; picomole detection levels: $\lambda_{\text{max}} = 254 \text{ nm}^{23\cdot24}$

<sup>\*</sup>EC = electrochemical; F = fluorescence; UV = ultraviolet; Vis = visible.

#### **Detection Reagents**

#### References

- 1. Durst, H.D., *et al.* (1975). *Anal. Chem.* **47,** 1797.
- 2. Borch, R.F., et al. (1975). Anal. Chem. 47, 2437.
- 3. Grushka, E. et al. (1975). J. Chromatogr. 112, 673.
- 4. Fitzpatrick, F.A. (1976). Anal. Chem. 48, 499.
- 5. Nagels, L., et al. (1980). J. Chromatogr. 190, 411.
- 6. Ahmed, M.S., et al. (1980). J. Chromatogr. 192, 387.
- 7. Pierce Technical Bulletin: Preparation of Phenacyl and p-Bromophenacyl Derivatives for HPLC.
- 8. Stocchi, V., et al. (1985). J. Chromatogr. 349, 77-82.
- 9. Chang, J.Y., et al.. (1981). Biochem. J. 199, 547-555.
- 10. Chang, J.Y., et al. (1982). Biochem. J. 203, 803-806.
- 11. Vendrell, J., et al. (1986). J. Chromatogr. 358, 401-413.
- 12. Chang, J.Y., et al. (1984). J. Chromatogr. 295, 193-200.
- 13. Lin J.K., et al. (1980). Anal. Chem. 52, 630-635.

- 14. Chang, J.Y., et al. (1981). FEBS Letters 132, 117-120.
- Marfey, P. (1984). Carlsberg Res. Comm. 49, 591-596.
   Jones, B.N. and Gilligan, J.P. (1983). Amer. Biotech. Lab. 12, 46-51.
- 17. Fiedler, H.P., et al. (1986). J. Chromatogr. 353, 201-206.
- 18. Fried, V.A., et al. (1985). Anal. Biochem. 146, 271-276.
- 19. Jones, B.N., et al. (1983). J. Chromatogr. 266, 471-482.
- 20. Seaver, S.S., et al. (Sept./Oct. 1984). Biotechniques 254-260.
- 21. Böhlen, P., et al. (1979). Anal. Biochem. 94, 313-321.
- 22. Stein, W.H. and Moore, S. (1949). Cold Spring Harbor Symp. Quant. Bio. 14, 179.
- 23. Heinrickson, R.L. and Meredith, S.C. (1984). Anal. Biochem. 137, 65-74.
- 24. Scholze, H. (1985). J. Chromatogr. 350, 453-460.
- 25. Caudill, W.L., et al. (1982). J. Chromatogr. 227, 331.
- 26. Caudill, W.L., et al. (1982). Bioanalytical System's Current Separations 4(4), 59.
- 27. Andrews, J.L., et al. (1982). Arch. Biochem. Biophys. 214, 386-396.
- 28. Aberhart, D.J., et al. (1985). Anal. Biochem. 151, 88-91.
- 29. Szokan, G., et al. (1988). J. Chromatogr. 444, 115-122.

#### FDAA, Marfey's Reagent

#### Makes it quick and easy for you to separate and quantitate optical isomers of amino acids by reverse-phase HPLC.

Optical isomers of amino acids can be simply and conveniently derivatized with Thermo Scientific FDAA, Marfev's Reagent (1-fluoro-2,4-dinitrophenyl-5-L-alanine amide) - and preparation is complete in just 90 minutes.

With Marfey's Reagent, the amino acid derivatives can easily be separated and quantitated by reverse-phase HPLC. Derivatives have an absorption coefficient of ~3 x 104 and can be detected by UV at 340 nm with picomole sensitivity.




Figure 2. Separation of D- and L-amino acids on 100 mm x 4.6 mm C18. Conditions: A) 0.05 M triethylamine phosphate, pH 3.0; B) acetonitrile. Linear gradient: 10 to 40% B in 45 minutes, 2.0 ml/minute, 25°C, 340 nm.

#### **PROTOCOL**

#### **Preparation of FDAA Derivatives**

- 1. Place 100  $\mu$ l (5  $\mu$ mol) sample in a 1.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Add 200 µl of a 1% (w/v) solution of FDAA in acetone. Add 40 µl of 1.0 M sodium bicarbonate. µmol FDAA: µmol amino acid should be
- 3. Heat at 40°C for 1 hour. Remove and cool.
- 4. Add 20  $\mu l$  2 M HCl. Allow sample to degas.
- 5. Analyze. Conditions: 100 mm x 4.6 mm C18; UV at 340 nm A: 0.05 M TEA phosphate, pH 3.0; B: CH3CN Linear gradient: 10% B to 40% in 45 minutes, Flow: 2.0 ml/minute at 25°C

1. Marfey, P. (1984). Carlsberg Res. Comm. 49, 591-596.

Ordering Information					
Product #	Description	Pkg. Size			
TS-48895	FDAA, Marfey's Reagent (1-fluoro-2,4-dinitrophenyl-5-L -alanine-amide)	50 mg			

#### p-Bromophenacylate Reagent

Procedure gives quantitative yields with few or no side reactions.

p-Bromophenacylate MW 277.94

Durst, et al. have described a novel preparation of various phenacyl esters and their use as UV visualizing agents in the 1-10 ng range. This procedure gives quantitative yields with few or no side reactions. Phenacyl esters have been used to separate many saturated and unsaturated fatty acids,23 including prostaglandins.4

Phenacyl esters have some significant advantages over previously reported methods, including:

- · Pre-mixing of phenacylbromide and crown ether is not necessarv
- Derivatization is both rapid and quantitative, with yields of more than 95% in 15-20 minutes at 80°C
- Excess reactants do not interfere
- Large excess of alkylating reagent is not necessary
- · Small amounts of water or alcohol do not interfere
- If isolation is desired, products usually are crystalline

#### **PROTOCOL**

#### **Preparation of Phenacyl Esters**

p-Bromophenacylate Reagent (0.1 µmol/ml p-Bromophenacylbromide, 0.005 µmol/ml crown ether in acetonitrile)

- 1. Dissolve ~10 mg acid in MeOH in a 5.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial fitted with Thermo Scientific Reacti-Vial Magnetic Stirrer. Neutralize to the phenolphthalein endpoint with KOH/MeOH.\*
- 2. Evaporate the MeOH with N<sub>2</sub>.
- 3. Add 1.0 ml Phenacylate Reagent and 2.0 ml dry CH<sub>3</sub>CN.
- 4. Heat at 80°C with stirring for 30 minutes.
- 5. Remove and cool.
- 6. Analyze. Conditions: C18; UV at 250 nm A: CH, CN; B: deionized H, 0 Linear gradient: 80% A to 100% A; Flow: 2.0 ml/minute
- \* If the formation of potassium salts is undesirable, neutralize by adding KHCO, at five times the total acid instead of using KOH.

Ordering Information					
Product #	Description	Pkg. Size			
TS-48891	p-Bromophenacylate Reagent 0.1 mmol/ml p-Bromophenacylbromide, 0.005 mmol/ml crown ether in acetonitrile	10 ml Hypo-Vial Sample Storage Vial			

# Ultimate core performance

- Rugged and reproducible solid core particles
- Fast separations with superb resolution
- Low backpressures





#### The History

The development of amino acid analysis began in 1820 when Braconnot isolated glycine from a hydrolyzate of gelatin.1 Later, in 1848, the Dutch chemist Mulder showed that glycine contains nitrogen, a major component of amino acids.<sup>2</sup> It was not until 1883, however, that Kjeldahl introduced a method that accurately determined the amount of nitrogen in a protein/amino acid sample.3

By 1910, most of the amino acids had been isolated and their structures discovered. As the number of known amino acids accumulated, it became possible to group them on the basis of common chemical features. At that time, it was discovered that all amino acids have the same general formula and differ only by the chemical structure of the side chains.

From 1910-1940, amino acid research was characterized by the work of quantitative analysts, as opposed to the organic chemists of the 1800s. Amino acid analyses conducted during the 1800s and early 1900s were laborious, often extending over weeks and months. While the amino acid content of a number of proteins was discovered, exact information was not always obtainable using the equipment available during that time period.4

The introduction of chromatography opened new doors for amino acid analysis. The first breakthrough came when Martin and Synge introduced partition chromatography, which separates the acetyl derivatives of certain amino acids.5 In this method, an equilibrium is established between two liquid phases. Silica gel is mixed with a solution of water and an indicator. The resulting slurry is packed into a column, forming the stationary phase. Next, the acetyl amino acids are dissolved in solvent, forming the mobile phase. The acids then are placed in the same column. The acetyl amino acids flow through the column at different rates. Separation is made visible by the bands of color change in the indicator.

While this system successfully separated mono-amino and mono-carboxylic acids, it was impractical for other types of amino acids.

Later, Martin and his associates used filter paper as an alternative to silica gel, developing a paper chromatography method that is still in use today. The amino acids were dissolved in butanol and allowed to seep onto the filter paper for a set amount of time. The paper then was dried and sprayed with a dilute solution of ninhydrin (2,2-dihydroxy-1,3-indandione) in butanol. The colored spots were measured and compared with the set values for those experimental conditions.

The separations achieved with paper chromatography were only semi-quantitative. Column chromatography, on the other hand, had potential for quantitation, but the separations were imperfect. The introduction of ion exchange chromatography solved these problems, allowing column separation of amino acids without any prior derivatization.<sup>6,7</sup> Initially used to remove carbohydrate contamination from starch columns, ion exchange resins were quickly found to have great potential for separating amino acids. While many types of polymeric exchange resins were tested, polysulfonic resins (such as Dowex® 50) provided the best separations.8

Modifications to these procedures have improved amino acid separations. Resin characteristics, column size, column temperature, buffer pH and ionic strength all have been modified to improve resolution of amino acid mixtures and achieve specific separations. Also, quantitation was greatly improved by the use of post-column reactions with ninhydrin. At one time ninhydrin was the most widely used detection system; however, more sensitive indicators, such as o-phthalaldehyde were developed to increase analytical sensitivity.

#### Developments in Amino Acid Analysis

Improvements in amino acid analysis by ion exchange chromatography have involved the analytical system, as well as the instrumentation. Systems have been developed (by varying buffer pH or ionic strength) that work to displace the amino acids into discrete bands. The buffer systems are compatible with single- or two-column analysis of amino acids found in protein hydrolyzates or physiological fluids. Buffer systems are determined by the counter ion used (sodium or lithium) and by the method of buffer changes introduced to the resin (step changes or gradient elution).9-15 The most commonly used buffering component, citrate, is suitable for solutions below pH 7.16 Buffers are prepared either with citric acid or an alkali salt, and citrate concentrations of 0.05 to 0.06 7 M are common.

Unfortunately, for high-sensitivity work, citric acid is a significant contributor to amino acid contamination. Therefore, to achieve consistent analyses, it is essential to use high-purity reagents for buffer preparation.

Alternatives to ion exchange are available for the separation of amino acids. Because amino acid analysis is one of the basic protein chemistry tools available, more rapid and sensitive methods for separation and quantitation are desirable.<sup>17</sup> Several pre-column derivatization methods using reverse-phase HPLC have been developed.

Two of the most widely used of these methods involve the formation of dansyl<sup>18-19</sup> or o-phthalaldehyde (OPA)<sup>20-24</sup> derivatives of amino acids prior to HPLC analysis. Both methods offer greater sensitivity and shorter analysis time than post-column derivatization techniques. Other methods include the quantitative derivatization of amino acids with phenylisothiocyanate (PITC) and the separation and quantitation of the resulting phenylthiocarbonyl derivatives via HPLC. These derivatives are stable enough to eliminate in-line derivatization.

# <u>Ultimate</u> **Bio HPLC performance**

- Industry leading columns
- High efficiency and resolution
- Excellent reproducibility

# ProPac, DNAPac, ProSwift columns



# **Thermo Scientific Web-Based Resources for Chromatographers**

Beyond this catalog, our Thermo Scientific Chromatography team shares its extensive expertise through our web-based Chromatography Resource Center and the Separated by Experience enewsletter.

The Chromatography Resource Center, accessed at www.thermoscientific.com/ chromatography, provides technical support, applications, technical tips and literature to help you move your separations forward, quickly and easily.

The bi-monthly Separated by Experience enewsletter keeps you up-to-date on the latest technical and product information of interest to chromatographers. Valuable information developed by chromatographers for chromatographers.

www.thermoscientific.com/chromatography

# **Separated by Experience**

As the world's sole manufacturer of Thermo Scientific Hypersil silica, we have set a very high standard in HPLC and continue to maintain it with innovative new products spanning across SPE, HPLC, GC and Vials & Closures. We offer one of the broadest selections of premier chromatographic phases and innovative hardware designs available, combined with superb technical support and customer service. Whether you use HyperSep SPE for rapid sample preparation, Hypersil GOLD for routine separations or are looking for something new for your most challenging methods like Accucore or SOLA or TraceGOLD GC columns, we have the choices to meet your needs.

# **TNBSA**

An excellent choice for spectrophotometric detection.

Trinitrobenzene sulfonic acid (TNBSA) reacts readily with the primary amino groups of amino acids in aqueous solution at pH8 to form yellow adducts. No colored derivatives are formed with the secondary amino acids proline and hydroxyproline. The colored derivatives are monitored at 345 nm and have extinction coefficients in the range of 1-1.5 x 104.

TNBSA has been used as a hydrophilic modifying reagent for the detection of primary amines in samples containing amino acids, peptides or proteins. It is an excellent reagent for rapid qualitative and quantitative estimation of these biomolecules.

# **Highlights:**

- · Couples with primary amines, sulfhydryls and hydrazides in aqueous solution at pH 8, without undesirable side reactions
- Excellent for solution or solid phase analysis
- Chromogenic ( $\lambda_{max} = 335 \text{ nm}$ )

## References

- 1. Goodwin, J.F. and Choi, S-Y. (1970) Clinical Chemistry. 16, 24-31.
- 2. Snyder, S.L. and Sobocinski, P.Z. (1975) Anal Biochem. 64, 284-288.
- 3. Drozdovskaya, N.R., et al. (1982) FEBS Lett. 150, 385
- 4. Takahashi, S., et al. (1984) Chem. Lett. (Jpn). 1, 127

Ordering	Information	
Product #	Description	Pkg. Size
<b>≭</b> TS-28997	TNBSA (2,4,6-Trinitrobenzene sulfonic acid; 5% w/v methanol solution)	100 ml

\* Additional hazardous handling charge.

Ninhydrin-based monitoring systems are among the most widely used methods for quantitatively determining amino acids after they are separated by ion exchange chromatography.

The color reaction between amino-containing compounds and ninhydrin (2,2-dihydroxy-1,3-indandione) is very sensitive. McCaldin has studied all phases of ninhydrin chemistry and proposed a mechanism for the reaction of ninhydrin with amino acids, accounting for the aldehydes, carbon dioxide, ammonia and hydrindantin known to be produced.<sup>1</sup> A yellow colored product (monitored at 440 nm) is formed upon reaction with the secondary amino acids, proline and hydroxyproline.2 Ninhydrin decarboxylates and deaminates the primary amino acids, forming the purple complex known as Ruhemann's Purple,3 which absorbs maximally at 570 nm.

Ninhydrin chemistry was adapted to a fully automatic, two-column amino acid analysis procedure in 1958 by Spackman, Stein and Moore. 4 Moore and Stein defined the requirements for a reducing agent (such as stannous chloride) to achieve reproducible color values for amino acids monitored with ninhydrin.<sup>5</sup> Titanous chloride was reported by James to eliminate precipitates encountered when using stannous chloride. 6-8 Methyl Cellosolve® (ethylene glycol monomethyl ether) buffered with 4 M sodium acetate at pH 5.51,9 and dimethylsulfoxide (DMSO) buffered with 4 M lithium acetate at pH 5.2010 are the most common solvents used for ninhydrin. DMSO remains stable longer than Methyl Cellosolve, particularly when kept chilled. These ninhydrin reagent solutions, with increased stability, were also reported by Kirschenbaum.11

Sensitivity of the ninhydrin system depends on several factors. Amino acids produce slightly different color yields, and these values may vary from one reagent preparation to the next. Ninhydrin also is sensitive to light, atmospheric oxygen and changes in pH and temperature. When ninhydrin becomes oxidized, its color does not develop well at 570 nm, but absorption at 440 nm remains fairly constant. When the height of the proline peak at 440 nm approaches the height of the glutamic acid peak at 570 nm, for equal amounts of each, reagent degradation is indicated.

#### References

- 1. McCaldin, D.J. (1960). Chem. Rev. 60, 39.
- 2. Hamilton, P.B. (1966). Advan. Chromatogr. 2, 3.
- 3. Ruhemann, J. (1911). J. Chem. Soc. London 99, 797.
- 4. Spackman, D.H., Stein, W.H. and Moore, S. (1958). Anal. Chem. 30, 1190.
- 5. Moore, S. and Stein, W.H. (1954). J. Biol. Chem. 211, 907.
- 6. James, L.B. (1971). J. Chromatogr. 59, 178.
- 7. James, L.B. (1978). J. Chromatogr. 152, 298-300.
- 8. James, L.B. (1984). J. Chromatogr. 284, 97-103.
- 9. Moore, S. and Stein, W.H. (1948). J. Biol. Chem. 176, 367.
- 10. Moore, S. (1968). J. Biol. Chem. 243, 6281-6283.
- 11. Kirschenbaum, D.M. (1965). Anal. Biochem. 12, 189.

Reaction Scheme. The course of the ninhydrin reaction with amino acids is as follows:

- 1. Ninhydrin (2,2-dihydroxy-1,3-indandione) reacted with amino acid.
- 2. The intermediate formed as the first reaction product.
- 3. Intermediate gives rise to dipolar ion by decarboxylation and dehydration.
- 4. The dipolar ion hydrolyzes, producing the amine.
- 5. The amine condenses with a second molecule of ninhydrin to give Ruhemann's Purple.

# Ninhydrin

The reagent of choice for detection of amino acids.

Since Stein and Moore pioneered amino acid chromatography in 1949,1 Thermo Scientific Ninhydrin has been used in amino acid chromatography advancements. The most recent techniques and sensitive instruments require the superb color response and a low blank that only our Ninhydrin offers.

Our Ninhydrin is indefinitely stable and requires no refrigeration or special care. Just keep the bottle tightly sealed, avoiding direct sunlight and ammonia in the laboratory atmosphere.

### References

- 1. Stein, W.H. and Moore, S. (1949). Cold Spring Harbor Symp. Quant Biol. 14, 179.
- 2. Moore, S. (1968). Amino acid analysis: aqueous dimethyl sulfoxide as solvent for the ninhydrin reaction. J. Biol. Chem. 243(23), 6281-6283.
- 3. James, L.B. (1978). Amino acid analysis: ninhydrin reaction with titanous chloride. J. Chromatogr. 152, 298-300.

Ordering	Information	
Product #	Description	Pkg. Size
TS-21003	Ninhydrin	500 g

# **National Scientific Mass Spec Certified Vials**

Industry's first and only pre-cleaned, low particle, low background chromatography vial

The only vial clean enough for every high sensitivity application

- Pre-cleaned vials provide unmatched consistency
- Pre-cleaned vial packaging protects the product integrity
- High purity closures packed in air-tight re-closeable container
- Tested and certified for 13 critical physical characteristics affecting vial performance

For your FREE sample pack, visit www.nationalscientific.com/cleanvials



The extraction and purification of proteins play an important role in determining amino acid content. These methods are based on one or more of their physical characteristics (e.g., solubility, molecular size, charge, polarity and specific covalent or noncovalent interactions). The techniques commonly used to separate proteins and peptides include:

- Reverse-phase HPLC
- · Polyacrylamide gel electrophoresis
- · Gel filtration
- Ion exchange chromatography
- · Affinity chromatography

Table 1 provides a more detailed list of methods for fractionating peptide mixtures.25



Table 1. Methods for the fractionation of peptide mixtures.

Technique	Properties of Peptide Molecules Exploited
Centrifugation	Solubility
Size exclusion chromatography	Size
Ion exchange chromatography	Charge, with some influence of polarity
Paper electrophoresis	Charge and size
Paper chromatography	Polarity
Thin layer electrophoresis	Charge and size
Thin layer chromatography	Polarity
Polyacrylamide gel electrophoresis	Charge and size
High-performance liquid chromatography (HPLC)	Polarity
Gas chromatography	Volatility of derivatives
Counter-current extraction	Polarity; sometimes specific interactions
Affinity chromatography	Specific interactions
Covalent chromatography or irreversible binding	Disulfide bond formation; reactivity of homoserine lactone

# Hydrolysis

Most protein samples require some form of chemical treatment before their component amino acids are suitable for analysis. Protein and peptide samples must be hydrolyzed to free amino acids from peptide linkages. Acids (usually HCI) are the most widely used agents for hydrolyzing proteins.

A simplified hydrolysis procedure involves refluxing the protein with excess HCI, then removing the excess acid in vacuum.26 The lyophilized protein then is suspended in constant boiling 6 N HCl and introduced into the hydrolysis tube. The sample is frozen by immersing the tube in dry ice and acetone. Before sealing, the tube is evacuated to avoid formation of cysteic acid, methionine sulfoxide and chlorotyrosine.27 This procedure minimizes decomposition of reduced S-carboxymethylcysteine and analyzes S-carboxymethylated proteins. Hydrolysis is conducted at 110°C (with the temperature accurately controlled) for 20-70 hours by Moore and Stein's method.<sup>28</sup> After hydrolysis, residual HCI is removed in a rotary evaporator. The residue is dissolved in water and brought to the appropriate pH for addition to the analyzer column.<sup>28</sup>

## References

- 1. Braconnot, H. (1820). Ann. Chem. Phys. 13, 113.
- 2. Mulder, G.J. (1848). J. Prakt. Chem. 38, 294.
- 3. Kjeldahl, J.G.C.T. (1883). Z. Anal. Chem. 22, 366.
- 4. Rattenbury, J.M. (1981). Amino Acid Analysis, Ellow Horwood, publisher, Chicester, England, Chapter 1.
- 5. Martin, A.J.P. and Synge, R.L.M. (1941). Biochem. J. 35, 1358.
- 6. Elsden, S.R. and Synge, R.L.M. (1944). Proc. Biochem. Soc. IX.
- 7. Synge, R.L.M. (1944). Biochem. J. 38, 285.
- 8. Stein, W.H. and Moore, S. (1949). Cold Spring Harbor Symposia Quant. Biol. 14, 179.
- 9. Liao, T.H., et al. (1973). Anal. Chem. 45, 2286.
- 10. Hare, P.E. (1972). Space Life Sciences 3, 354
- 11. Dus, K., et al. (1966). Anal. Biochem. 14, 41.
- 12. Thomson, A.R. and Miles, B.J. (1964). Nature 203, 483.
- 13. Robinson, G.W. (1975). J. Chromatogr. 3, 416.
- 14. Benson, J.V. (1972). Anal. Biochem. 50, 477.
- 15. Benson, J.V. and Patterson, J.A. (1965). Anal. Chem. 37, 1108.
- 16. Bates, R.G. and Pinching, G.D. (1949). *J. Amer. Chem. Soc.* **71**, 1274. 17. Jones, B.N. and Gilligan, J.P. (1983). *American Biotech. Lab.* December, 46-51.
- 18. Weidmeier, V.T., et al. (1982). J. Chromatogr. 231, 410-417 19. DeJong, C., et al. (1982). J. Chromatogr. 241, 345-359.
- 20. Lindroth, P. and Mopper, K. (1979). *Anal. Chem.* **51**, 1667-1674. 21. Jones, B.N., *et al.* (1981). *J. Liq. Chromatogr.* **4**, 565-586.
- 22. Turnell, D.C. and Cooper, J.D.H. (1982). Clin. Chem. 28, 527-531.
- 23. Umagat, H., et al. (1982). J. Chromatogr. 241, 324-334.
- 24. Jones, B.N. and Gilligan, J.P. (1983). *J. Chromatogr.* **266**, 471-482.
- 25. Allen, G. (1981). Laboratory Techniques in Biochemistry and Molecular Biology, T.S. Work and R.H. Burdon, Eds., Elsevier/North-Holland, Biomedical Press.
- 26. Moore, S. and Stein, W.H. (1963), Meth. Enzymol. 6, 819.
- 27. Eveleigh, J.W. and Winter, G.D. (1970). Protein Sequence Determination, S.B. Needleman, Ed., Springer-Verlag, pp. 92-95.
- 28. Moore and Stein, op. cit.

# Constant Boiling (6N) Hydrochloric Acid

## For total protein hydrolysis.

Thermo Scientific Hydrochloric Acid is purified and packaged to ensure a ninhydrin negative blank on hydrolysis. Convenient, pre-scored ampule packaging of the ready-touse HCI maintains reagent integrity. This virtually eliminates exposure to laboratory atmospheres, fingerprints and other contaminants resulting from pipetting from bulk bottles.

Eveleigh and Winter give an excellent description of the total protein hydrolysis technique using Constant Boiling Hydrochloric Acid.1 Standard protein hydrolysis conditions are 105-110°C for 16-24 hours. At 150°C, this reagent can hydrolyze peptides in 6 hours.

# **Highlights:**

- Hydrolyzes peptides in 6 hours at 150°C
- Specially purified to give ninhydrin-negative blank on hydrolysis
- Packaged in ampules to eliminate contamination and ensure product integrity

#### References

- 1. Eveleigh, J.W. and Winter, G.D. (1970). Protein Sequence Determination, Ed Needleman, S.B., Springer-Verlag, pp 92-95.
- 2. Blankenship, et al. (1989). High-sensitivity amino acid analysis by derivatization with o-Phthaldialdehyde and 9-Fluorescence detection: applications in protein structure determination. Anal. Biochem. 178, 227-232.
- 3. Hurley, J.B., et al. (1984). Isolation and characterization of a cDNA clone for the subunit of bovine retinal transducin. Proc. Natl. Acad. Sci. U.S.A. 81, 6948-6952.
- 4. Lee, K., et al. (1979). Derivatization of cysteine and cystine for fluorescence amino acid analysis with the o-Phthaldialdehyde/2-mercaptoethanol reagent. J. Biol. Chem. July 25, 6248-6251.

## **Ordering Information**

Product #	Description	Pkg. Size
TS-24308	Hydrochloric Acid [Constant boiling, Hydrochloric Acid 6N Sequencing Grade]	10 x 1 ml ampules

# Amino Acid Standard H

# Our high-purity amino acid calibration standard for protein hydrolyzates.

The high-purity amino acids of Thermo Scientific Amino Acid Standard H are ideal for calibrating amino acid analyzers. To permit standardization of microbiological and other assays, we have used the L-form configuration. The molar concentration of each standard is verified by conventional amino acid analysis methods.

With the exception of cystine, each amino acid is supplied at a concentration of 2.5 µmoles/ml in 0.1 N HCI. The following amino acids are included with our Amino Acid Standard H:

<ul> <li>L-Alanine</li> </ul>	L-Leucine
<ul> <li>Ammonia [(NH<sub>4</sub>)<sub>2</sub> SO<sub>4</sub>]</li> </ul>	L-Lysine•HCI
• L-Arginine	L-Methionine
<ul> <li>L-Aspartic Acid</li> </ul>	L-Phenylalanine
<ul> <li>L-Cystine</li> </ul>	L-Proline
<ul> <li>L-Glutamic Acid</li> </ul>	L-Serine
• Glycine	L-Threonine
<ul> <li>L-Histidine</li> </ul>	L-Tyrosine
L-Isoleucine	L-Valine

## **Instructions for Use**

Thaw Standard H and shake well. Dilute appropriately with suitable buffers to a concentration compatible with the full-scale sensitivity of your amino acid analyzer.

## **Storage**

When kept frozen, an unopened vial has an indefinite storage life. Once the seal is broken, the reagent has a maximum storage life of six months. Store Amino Acid Standard H frozen between uses.

Ordering	<b>Information</b>	
Product #	Description	Pkg. Size
TS-20088	Amino Acid Standard H	10 x 1 ml ampules

# PITC (Phenylisothiocyanate)

Ideal for the quantitative pre-column derivatization of amino acids by reverse-phase HPLC.14

Thermo Scientific PITC, also known as Edman's Reagent, reacts readily with amino acids in 5-10 minutes at room temperature. The resulting phenylthiocarbamyl derivatives can be separated and quantified in 30 minutes using reverse-phase HPLC. This method produces stable products with all amino acids, including proline.

# TO COUPLE AMINO ACID STANDARD H WITH PITC.1

- 1. Dry 10 µl Amino Acid Standard H in a small test tube. Dissolve dried standard in 100 µl coupling buffer (acetonitrile: pyridine: triethylamine:
- 2. Dry standard solution by rotary evaporation. Dissolve the residual amino acids once more in 100 µl of coupling buffer.
- 3. Add 5 µl of PITC.
- 4. Allow a 5-minute reaction at room temperature.
- 5. Evaporate sample to dryness by rotary evaporation under high vacuum.
- 6. Dissolve the resulting PITC-amino acids in 250 µl of 0.05 M ammonium acetate, water or water:acetonitrile (7:2).
- 7. Analyze quantities of 1 to 10  $\mu$ l (100 to 1,000 picomoles of each amino acid) by reverse-phase HPLC.

NOTE: Make certain that all HCl is evaporated before derivatization.

# References

- 1. Heinrikson, R.L. and Meridith, S.C. (1984). Anal. Biochem. 136, 65-74.
- 2. Scholze, H. (1985). J. Chromatogr. 350, 453-460.
- 3. Janssen, et al. (1986). Chromatogr. 22(7-12).
- 4. Evert, R.F. (1986). Anal. Biochem. 154, 431-435.

Ordering	g Information	
Product #	Description	Pkg. Size
TS-26922	PITC (Phenylisothiocyanate)	10 x 1 ml ampules
TS-20088	Amino Acid Standard H	10 x 1 ml ampules

# **Dabsyl Chloride**

# It's recrystallized twice for twice the quality!

$$\begin{array}{c|c} & & & & & & & & & & & & & & & & \\ & & & & & & & & & & & & & \\ & & & & & & & & & & \\ & & & & & & & & & \\ & & & & & & & & \\ & & & & & & & & \\ & & & & & & & \\ & & & & & & & \\ & & & & & & & \\ & & & & & & \\ & & & & & & \\ & & & & & & \\ & & & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & & & & \\ & & & & \\ & & & & & \\ & & & & \\ & & & & \\ & & & & \\ & & & \\ & & & \\ & & & \\ & & & \\ & & & \\ & \\$$

Dabsyl Chloride MW 323.80

Thermo Scientific Dabsyl Chloride is ideal for the pre-column derivatization and detection of amino acids in visible light down to sub-picomolar levels, followed by reverse-phase HPLC.

# **Highlights:**

- Analysis of 10-30 ng of protein hydrolyzates<sup>1,2</sup>
- Analysis of peptides and determination of C-terminal sequence of polypeptides<sup>1</sup>
- Analysis of phospho-amino and amino acid amides<sup>3</sup>
- Analysis of amino acid neurotransmitters in mouse brain4
- Optimal reaction conditions<sup>5</sup>

#### References

- 1. Chang, J.Y., et al. (1981). Biochem. J. 199, 547-555.
- 2. Chang, J.Y., et al. (1982). Biochem. J. 199, 803-806
- 3. Chang, J.Y. (1984). J. Chromatogr. 295, 193-200.
- 4. Chang, J.Y., et al. (1981). FEBS Lett. 132, 117-120.
- 5. Vendrell, J., et al. (1986). J. Chromatogr. 358, 401-413.
- 6. Lin, J.K., et al. (1980). Clin. Chem. 26, 579-583
- 7. Chang, J.Y., et al. (1983). Methods. Enzymol. 92, 41-48.
- 8. Stocchi, V., et al. (1985). J. Chromatogr. 349, 77-82.

Ordering Information		
Product #	Description	Pkg. Size
<b>≭</b> TS-21720	Dabsyl Chloride (4-Dimethylaminoazobenzene -4-Sulfonyl Chloride)	500 mg

\* Additional hazardous handling charge.

# FDAA, Marfey's Reagent

Derivatizes optical isomers of amino acids in just 90 minutes.

Thermo Scientific FDAA, (1-fluoro-2,4-dinitrophenyl-5-Lalanine amide), offers complete derivatization of amino acid isomers in 90 minutes. Derivatized amino acids then are separated and quantitated by reverse-phase HPLC. The nature of the reagent and the resultant reaction products with D-diastereomers suggest that strong intramolecular hydrogen bonding causes these derivatives to elute much later than their L-diastereomer counterparts. Derivatives have an absorption coefficient of approximately 3 x 104. They can be detected by UV at 340 nm with picomole sensitivity. Complete instructions are included with each order.

### TO PREPARE FDAA DERIVATIVES:

- 1. Place 100 µl (5 µmoles) sample in a 1.0 ml Thermo Scientific Reacti-Vial Small Reaction Vial.
- 2. Add 200 µl of 1% (w/v) solution of FDAA in acetone. Add 40 µl of 1.0 M sodium bicarbonate (µmoles FDAA: µmoles amino acid should be 1.5:1.0.)
- 3. Heat at 40°C for 1 hour. Remove and cool.
- 4. Add 20 µl 2 M HCI. Allow sample to degas.
- 5. Analyze.

Conditions: Hypersil GOLD 100mm x 4.6mm (part number 25005-104630) UV at 340 nm A: 0.05 M TEA phosphate pH 3.0 B: CH CN

Linear Gradient, 10% B to 40% B in 45 minutes Flow: 2.0 ml/min. at 25°C

- 1. Marfey, P., et al. (1984). Carlsberg Res. Comm. 49, 585-590.
- 2. Marfey, P. (1984). Carlsberg Res. Comm. 49, 591-596.
- 3. Szókán, G., et al. (1988). Applications of Marfey's Reagent in racemization studies of amino acids and peptides. J. Chromatogr. 444,115-122.
- 4. Aberhart, D.J., et al. (1985). Separation by high-performance liquid chromatography of (3R)- and (3S)-B-Leucine as diastereomeric derivatives. 1151, 88-91.
- 5. Martinez del Pozo, et al. (1989). Stereospecificity of reactions catalyzed by bacterial D-amino acid transaminase. J. Biol. Chem. 264(30) 17784-17789.

Ordering	J Information	
Product #	Description	Pkg. Size
TS-48895	FDAA Marfey's Reagent (1-fluoro-2, 4-dinitrophenyl-5-L -alanine amide)	50 mg

# Ultra-Pure Solvents for Amino Acid **Analysis**

# Ideal for HPLC and spectrophotometric applications.

Thermo Scientific HPLC Grade Solvents and Water are specially purified by proprietary methods and tested to ensure lot-to-lot consistency with a low UV absorbance to provide you with the most sensitive detection across all wavelengths. All are packaged in amber glass bottles and sealed with Teflon TFE-lined fluorocarbon caps for ultimate protection. Our solvents are then tested to the highest specifications to ensure the integrity of your data, maximized sensitivity in your assay and prolonged life of your equipment.

# **Physical Properties Acetonitrile. HPLC Grade**

• UV Cutoff: 190 nm

• Optical Absorbance: <0.02 at 220 nm

• Refractive Index at 25°C: 1.342

# Water, HPLC Grade

• UV Cutoff: 190 nm

• Optical Absorbance: <0.005 at 220 nm

Refractive Index at 25°C: 1.332

# Dimethylformamide (DMF), Sequencing Grade

• HCON (CH<sub>2</sub>)<sub>2</sub>

• Purity (GC): ≥99%

• MW: 73.09

• Density: 0.944

• B.P. 153°C

Water: <0.1%</li>

## Dimethylsulfoxide (DMSO), Sequencing Grade

• C.H.OS

• Purity (GC): >99.5%

• MW: 78.13 • Density: 1.101 • Water: ≤0.2%

## **Pyridine**

• C<sub>c</sub>H<sub>c</sub>N

• Purity (GC): ≥99%

• MW: 79.10

• Density: 0.978

• B.P.: 115°C

Product #	Description	Pkg. Size
<b>♦ TS-51101</b>	Acetonitrile	1 L
TS-51140	Water	1 L
TS-20673	Dimethylformamide (DMF)	50 ml
<b>♦ TS-25104</b>	Pyridine (C <sub>E</sub> H <sub>E</sub> N)	100 g
<b>≮</b> TS-20688	Dimethylsulfoxide (DMSO)	950 ml

\* Additional hazardous handling charge. For GC Grade DMF and DMSO, see page 14.

# **Aspire Protein A** & Protein G Tips

# **Purify antibodies in minutes**

- Capture and purify a wide range of monoclonal and polyclonal IgG antibodies
- Embedded with high quality immobilized Thermo Scientific Protein A and Protein G Plus Agarose with the capacity to purify (1mg of human IgG)
- Spin column capacity for a fraction of the cost
- Fast and easy protocol does not compromise purity and yield



# Peptide Retention Time Calibration Mixture

# Heavy peptide mixture for column assessment and prediction of peptide retention times.

- Assessment of chromatography and MS instrument performance
- Prediction of peptide retention across multiple instrument platforms
- · Prediction of peptide retention time from sequence using calculated hydrophobicity factor
- Optimization of scheduled MS acquisition windows for improving quantification and increased multiplexing
- Internal standard to normalise for variation in retention times and peak intensities between runs

Ordering Information		
Product #	Description	Pkg. Size
TS-88320	Peptide Retention Time Calibration Mixture, 0.5pmol/µL	50μL x1
TS-88321	Peptide Retention Time Calibration Mixture, 5pmol/µL	200μL x1

# Peptide Retention Standard for Reverse-Phase HPLC

# Increases the efficiency of peptide elution profile predictions.

A simple, quantitative method for predicting peptide retention times was developed by Guo, et. al. 1-3 Retention times are predicted by totaling the values that represent the contribution in minutes of each amino acid residue and the peptide terminal groups.

Retention time is dependent upon the molecular weight of the peptide. The effect on retention is relatively unimportant with a small peptide, but it increases with the size of the molecule. The accuracy of predicting peptide retention time significantly decreases beyond 20 residues.

To ensure accuracy, a peptide standard is used to correct for instrument variation, column aging, n-alkyl chain length variation and ligand density.

By using a Thermo Scientific Peptide Retention Standard,

- Determine the relative order of peptide elution of a complex mixture
- · Increase the efficiency of predicting peptide elution profiles
- Save time in peptide purification
- Simplify identification of specific peptides in a complex
- Predict the HPLC retention time for peptides of known amino acid composition on reverse-phase HPLC columns
- Monitor column performance efficiency, selectivity and resolution during column aging
- · Compare reverse-phase columns from different manufacturers
- Evaluate reverse-phase supports of varying n-alkyl chain lengths and ligand densities

## References

- 1. Guo, D., et al. (1985). Proceedings of the Ninth American Peptide Symposium, Published by Pierce, Rockford, Illinois, page 23.
- 2. Guo, D., et al. (1986). J. Chromatogr. 359, 499-517.
- 3. Guo, D., et al. (1986). J. Chromatogr. 359, 519-532.
- 4. Mant, C.T. and Hodges, R.S. (1986). L.C. Magazine Liq. Chrom. and HPLC 4(3),
- 5. Guo, D., et al. (1987). J. Chromatogr. 386, 205-222.

Ordering	J Information	
Product #	Description	Pkg. Size
TS-31700	Peptide Retention Standard, S1-S5 Contains: 5 C-terminal amide decapeption	1 vial
	4 of which are N-acetylated w sequence variation as follows AC-Arg-Gly-X-X-Gly-Leu-Gly-I	vith the ::
	Amide, Gly³-Gly⁴-Ala³-Gly⁴, Ýa Val⁴. The fifth peptide, Ala³-Gl∙ N-amino group. This mixture v 100-200 injections at 0.1 AUFS	y⁴, contains a free vill provide



Note: Reacti-Therm Units are NOT supplied with a Reacti-Vap and Reacti-Blocks. These items must be purchased separately

# Reacti-Therm™ Dry Block Sample Incubation System

Featuring our exclusive Reacti-Therm Modules for dry block heating - and no sample contamination!

The Thermo Scientific Reacti-Therm System delivers uniform dry heat with unmatched convenience and versatility. The dry heat prevents many of the problems associated with water baths, including sample contamination.

Reacti-Therm Modules are easy-to-use, constant-temperature heaters that are ideal for your routine incubations. They also provide constant temperature control for samples held in Thermo Scientific Reacti-Vial Small Reaction Vials, as well as for samples in test tubes, microcentrifuge tubes and other small containers. Most applications that require heating, stirring or evaporation of small samples would benefit from the convenience and efficiency of Reacti-Therm Modules. These applications include:

- · Sample incubation
- Sample evaporation
- · Protein hydrolysis
- Small-scale reactions
- · Vacuum hydrolysis for amino acid analysis
- Derivatization reactions for HPLC and GC

Our Reacti-Therm Modules transfer heat through an aluminum alloy block. They hold a wide variety of interchangeable Thermo Scientific Reacti-Block Aluminum Blocks (page 46-47). Choose from four module designs to meet your exact incubation needs.

Reacti-Therm Heating Modules and Reacti-Therm Heating/ Stirring Modules – now quicker than ever!

Reacti-Therm Heating Modules feature a solid state electronic control. This highly efficient control system allows faster and easier temperature settings.


Dual Volt Modules				
Product #	Description	Pkg. Size		
TS-18820	Reacti-Probe Remote Temperature Probe	1 unit		
TS-18821	Reacti-Therm Heating/Stirring Module (Single Block)	1 unit		
TS-18822	Reacti-Therm Heating Module (Single Block)	1 unit		
TS-18823	Reacti-Therm Heating/Stirring Module (Triple Block)	1 unit		
TS-18824	Reacti-Therm III Heating Module (Triple Block)	1 unit		
TS-18825	Reacti-Vap Evaporator	1 unit		
TS-18826	Reacti-Vap Evaporator	1 unit		

Underwriters Laboratories, Inc. Listed

Note: Our Reacti-Therm Modules bear a CE marking for meeting the requirements of the European Union's Low-Voltage and EMC Directives.

# Reacti-Vap Evaporators

# Sample evaporation made easy!





Thermo Scientific Reacti-Vap Evaporator (9-port) and Reacti-Vap III Evaporator (27-port).

Thermo Scientific Reacti-Vap Evaporators are precisionmachined gassing manifolds. They provide simple, efficient evaporation by allowing the simultaneous or separate delivery of nonreactive pressurized gas to samples.

The Reacti-Vap III Evaporator triples the number of samples you can evaporate. Nine needles attach to each of the three individually regulated chambers. The evaporating head tilts back for easy needle attachment and removal.

The standard Reacti-Vap Evaporator attaches easily to single-block Reacti-Therm Modules. The Reacti-Vap III unit attaches easily to Reacti-Therm III Modules.

Ordering	J Information
Product #	Description
TS-18825	Reacti-Vap Evaporator (9-port) For use with Reacti-Therm Single Block Modules; TS-18822 and TS-18821, Includes 9 needles and plugs
TS-18826	Reacti-Vap III Evaporator (27-port) For use with Reacti-Therm III Modules; TS-18823 and TS-18824, Includes 27 needles and plugs

# Reacti-Vap Standard and Teflon Coated Needles

# Reduce cross-contamination and corrosion.



Thermo Scientific Reacti-Vap Teflon Coated Needles are made exclusively for use in Reacti-Vap Evaporators. They are blunt-ended, 19-gauge, stainless steel needles that reduce cross-contamination and corrosion when evaporating solvents that contain strong acids.

Each Reacti-Vap Needle has a Luer-Lok® hub for leak-proof attachment to Reacti-Vap Evaporators Needles are available in 4- and 6-inch lengths.

Ordering	Ordering Information					
Product #	Description	Pkg. Size				
TS-18782	Reacti-Vap Replacement Tube Kit 2.5 inch (64 mm)	Pkg. of 9 and plugs				
TS-18784	Reacti-Vap Teflon Coated Needles 4-inch (102 mm) x 19 gauge	Pkg. of 9				
TS-18786	Reacti-Vap Teflon Coated Needles 6-inch (152 mm) x 19 gauge	Pkg. of 9				

# Reacti-Vap Replacement Parts

Ordering Information					
Product #	Description	Pkg. Size			
TS-18827	Replacement Luer-Lok Fitting	Pkg. of 1			
TS-18828	Replacement Screws for Mounting Bracket	Pkg. of 4			
TS-18829	Replacement Height Adjustment Knob	Pkg. of 1			
TS-18830	Replacement Mounting Bracket	Pkg. of 1			
TS-18831	Replacement Metal Rod	Pkg. of 1			
TS-18832	Replacement Dial for Flow Control	Pkg. of 1			
TS-18833	Replacement Long Screws for Mounting Bracket	Pkg. of 4			

# Reacti-Therm Thermometers

## Teflon-coated, designed specifically for dry incubations.

Ordering	Information
Product #	Description
TS-18914	Reacti-Therm Thermometer, Mercury-free (0-100°C)
TS-18915	Reacti-Therm Thermometer, Mercury-free (0-200°C)

# Reacti-Block Aluminum Blocks

# There is one that is right for your sample needs!



Thermo Scientific Reacti-Block Aluminum Blocks are available with many hole configurations, machine-drilled to accommodate almost any size Reacti-Vial Small Reaction Vial (page 48), test tube or microcentrifuge tube. These highly efficient

units are constructed of an aluminum alloy for optimal thermal conductivity. To ensure proper heat transference, be sure to have a close block-to-sample container fit.

Each Reacti-Block Aluminum Block contains a thermometer well 7.1 mm dia. x 36.5 mm deep (excluding blank block J and K). Block dimensions are 9.4 cm long x 7.5 cm wide x 5.1 cm tall for all blocks except F, G, J and M which are 9.4 cm long x 7.5 cm wide x 7.6 cm tall.

The following Reacti-Block Aluminum Blocks can be used with all Reacti-Therm Modules including those equipped with Reacti-Vap Evaporators. Blocks B-1 and T-1 are specifically designed for use with Reacti-Vap Units.

# To complete your Reacti-Therm System Order:

- 1. Reacti-Therm Module
- 2. Reacti-Block Aluminum Block(s)
- 3. Reacti-Therm Thermometer
- 4. Reacti-Vap Evaporator

Ordering In	formation	
Product #		Description
TS-18801		Reacti-Block A-1 Holds 13 x 0.3 ml or 1 ml Reacti-Vials™; 13 holes/14 mm dia. x 23 mm deep
TS-18802		Reacti-Block B-1 Holds 9 x 3 ml or 5 ml Reacti-Vials; 9 holes/21 mm dia. x 32 mm deep
TS-18803		Reacti-Block C-1 Holds 13 x 3.5 ml Screw Cap Septum Vials; 13 holes/15 mm dia. x 34 mm deep
TS-18804		Reacti-Block Z-1 Holds 9 x 0.1mm Reacti-Vials; 9 holes/12mm dia. x 21mm deep
TS-18811		Reacti-Block M-1 Holds 6 x 27.5 ml Reacti-Vials; 6 holes/28.5 mm dia. x 70 mm deep
TS-18816		Reacti-Block S-1 Holds 13 x 13 mm dia. Test Tubes; 13 holes/14 mm dia. x 45 mm deep
TS-18817		Reacti-Block T-1 Holds 9 x 16 mm dia. Test Tubes; 9 holes/17 mm dia. x 45 mm deep
TS-18818		Reacti-Block U-1 Holds 8 x 20 mm dia. Test Tubes; 8 holes/21 mm dia. x 45 mm deep
TS-18819		Reacti-Block V-1 Holds 17 Microcentrifuge Test Tubes; 17 holes/11 mm dia. x 45 mm deep

# Reacti-Block Aluminum Blocks (continued)

The Reacti-Block Aluminum Blocks featured below are designed to be used exclusively with the Reacti-Therm Modules. The hole patterns do not match the needle configuration of Reacti-Vap Evaporators.

Ordering Int	formation		
Product #		Description	
TS-18806		<b>Reacti-Block F</b> Holds 8 x 6 ml Vacuum Hydrolysis Tubes; 8 holes/10 mm dia. x 64 mm deep	
TS-18807	• •	<b>Reacti-Block G</b> Holds 4 x 18 ml Vacuum Hydrolysis Tubes; 4 holes/19 mm dia. x 64 mm deep	
TS-18809		Reacti-Block J Blank/no holes (for custom drilling) 7.6 cm tall	
TS-18810		Reacti-Block K Blank/no holes (for custom drilling) 5.1 cm tall	
TS-18812		Reacti-Block L Holds 16 x 0.1 ml Reacti-Vials; 16 holes/12 mm dia. x 21 mm deep	

# Reacti-Vial Magnetic Stirrers

Offer faster reaction times with smooth mixing of small samples.



Mounted on a triangular matrix, these small Teflon-coated stirring bars fit the cone portion of 0.3, 1.0, 3.0, 5.0 and 10.0 ml Thermo Scientific Reacti-Vial Small Reaction Vials. For more information on Reacti-Vial Small Reaction Vials, see page 48.

Product #	Description	Pkg. Size
TS-16000	Reacti-Vial Magnetic Stirrers For use with 3.0, 5.0 and 10 ml Reacti-Vial Small Reaction Vials	Pkg. of 6
TS-16010	Reacti-Vial Magnetic Stirrers For use with 0.3 and 1.0 ml Reacti-Vial Small Reaction Vials	Pkg. of 6

When used with Thermo Scientific Reacti-Therm Heating/ Stirring Modules, these efficient stirrers provide:

- Faster reaction times with smooth, efficient mixing of small reaction samples
- Solubilization of sticky concentrated residues such as those found on evaporation of sugar solutions
- Increased speed-of-surface reactions by keeping insoluble reactants in suspension

# Reacti-Vial Small Reaction Vials

# Make small-sample handling easy and convenient.



# **Ideal for:**

- Residue isolation
- Derivative preparation
- Maximum sample retrieval
- Moisture protection
- Sample storage

Thermo Scientific Reacti-Vial Small Reaction Vials have an internal cone designed to make small-sample collection and handling easy and convenient. The cone feature is particularly useful for removing small quantities of sample with a syringe, even into the microliter range. The extra thick glass wall magnifies the sample, making these units ideal for observing chemical reactions. Reacti-Vial Small Reaction Vials can be used for derivatization, isolation and purification. You can also use Reacti-Vial Small Reaction Vials for precipitations, centrifugations and solvent separations.

Our amber Reacti-Vial Small Reaction Vials are manufactured from amber glass, and are amber throughout. These amber vials assure that your light-sensitive compounds are well protected. All Reacti-Vial Small Reaction Vials are supplied complete with Open-Top Screw Caps and Teflon/Rubber Laminated Discs (other discs can be ordered separately, see optional accessories below).

## **Reacti-Vial Small Reaction Vials**

			Inside Diameter Thread (mm) Style	Clear Pkg. of 12	Amber Pkg. of 12  Product #
Size	Dimensions (Diam. x Height) (mm ± 1mm)	Diameter		Product #	
100 μΙ	12 x 32	8	425-8	TS-13100	
0.3 ml	13 x 32	11	425-13	TS-13220	
1.0 ml	13 x 45	11	425-13	TS-13221	TS-13097
3.0 ml	20 x 47	18	425-20	TS-13222	
5.0 ml	20 x 60	18	425-20	TS-13223	TS-13099
10.0 ml	25 x 69	22	425-24	TS-13225	

# **Optional Accessories**

	Teflon/ Silicone Discs Pkg. of 72	Rubber Laminated Discs Pkg. of 72	Open-Top Screw Caps Pkg. of 72	Miniert Valves Pkg. of 72	Reacti-Vial Magnetic Stirrers Pkg. of 6
Vial Size	Product #	Product #	Product #	Product #	Product #
100 μΙ	TS-12708		TS-13208		
0.3 ml	TS-12712	TS-12412	TS-13215		TS-16010
1.0 ml	TS-12712	TS-12412	TS-13215		TS-16010
3.0 ml	TS-12718	TS-12418	TS-13218	TS-10135	TS-16000
5.0 ml	TS-12718	TS-12418	TS-13218	TS-10135	TS-16000
10.0 ml	TS-12722	TS-12422	TS-13219	TS-10130	TS-16000

<sup>\*</sup>For more information about discs and Mininert™ Valves, see page 50, 51.

<sup>\*</sup>For more information about Reacti-Vial Magnetic Stirrers, see page 47.

# Screw Cap Septum Vials

Autoclavable, borosilicate glass available in clear or amber.



- Sample collection the 40 ml clear or amber vial with a Thermo Scientific Disc is suitable for discrete water sampling under EPA 40 CFR Parts 136 and 141
- Automated GCs and LCs the 1.5 ml clear and amber vials fit autosamplers using standard 12 x 32 mm vials
- Heavy-duty, flip-top divider box provides easy access to vials; caps and septa and offers a convenient sample storage center

Thermo Scientific Screw Cap Septum Vials are supplied complete with Open-Top Screw Caps.

For economy, convenience and versatility in a vial and closure system, our Screw Cap Septum Vials are your best choice. A wide assortment of special closures and accessories make this system perfect for:

- Storage of reagents and standards under complete seal with instant syringe access
- Small derivatization reactions

## **Screw Cap Septum Vials**

				<b>Clear</b> Pkg. of 12	Amber Pkg. of 12
Size	Dimensions (Diam. x Height) (mm ± 1mm)	Inside Diameter (mm)	Thread Style	Product #	Product #
1.5 ml	12 x 32	8	425-8		TS-13080
3.5 ml	15 x 45	12	425-13	TS-13019	
7.0 ml	17 x 60	13	425-15	TS-13028	
14.0 ml	21 x 70	16	425-18	TS-13043	

## **Optional Accessories**

	Teflon/ Silicone Discs Pkg. of 72	Teflon/ Rubber Laminated Discs Pkg. of 72	Mininert Valves Pkg. of 72	Open-Top Screw Caps Pkg. of 72
Vial Size	Product #	Product #	Product #	Product #
1.5 ml	TS-12708			TS-13208
3.5 ml	TS-12712	TS-12412		TS-13215
7.0 ml	TS-12713			TS-13216
14.0 ml	TS-12716			TS-13217
40.0 ml	TS-12722	TS-12422	TS-10130	TS-13218

For more information about discs and Mininert Valves, see page 50, 51.

# Teflon/Silicone Discs

Unique discs that combine the inertness of a Teflon coating with the resealability of silicone.



Autoclavable Thermo Scientific Teflon/Silicone Discs are specifically designed to combine the resealability of silicone with the inertness of a Teflon coating. Many sizes are available to fit our Hypo-Vial and Reacti-Vial Screw Cap.

Structurally bonded (not cemented) Teflon coating to silicone. No cement to be leached or baked out of your sample after needle penetration.

# **Highlights:**

- · Reseals instantly puncture after puncture
- Compresses, giving it a "lock-washer" effect in maintaining a tight seal, and forces the Teflon coating to conform to the sealing surface
- No bent needles from a septum that is too hard; standard syringe and CC needles penetrate the entire disc with
- Protection of Teflon coating: The Teflon layer is a full 10 mils thick

# Thermo Scientific Teflon/Silicone Discs (pkg. of 72)

Product #	Silicone Thickness (mils)	Teflon Thickness (mils)	Fits These Containers
TS-12708	75	10	100 µl Reacti-Vial Small Reaction Vials 1.5 ml; Screw Cap Septum Vials
TS-12712	90	10	0.3 and 1.0 ml Reacti-Vial Small Reaction Vials 3.5 ml Screw Cap Septum Vials
TS-12713	90	10	7 ml Screw Cap Septum Vials
TS-12716	90	10	14 ml Screw Cap Septum Vials
TS-12718	90	10	3 and 5 ml Reacti-Vial Small Reaction Vials
TS-12720	125	10	6-125 ml Hypo-Vial Sample Storage Vials
TS-12722	90	10	25 and 40 ml Screw Cap Septum Vials 10 ml Reacti-Vial Small Reaction Vials

# Teflon/Rubber Laminated Discs

For a highly inert and nonreactive seal.



Thermo Scientific Teflon/Rubber Laminated Discs are constructed of white pharmaceutical rubber with 5 mils of Teflon coating bonded to one side. Total thickness of the disc is approximately 60 mils. The discs are excellent for use as cap liners when a highly inert nonreactive seal is desired.

Our rigid Teflon/Rubber Laminated Discs are more difficult to puncture than Teflon/Silicone Discs. Consequently, care must be taken when puncturing these discs to avoid bending the needle.

Our Teflon/Rubber Laminated Discs are auto-clavable and demonstrate no loss of integrity after heating above 100°C for 5 hours.

## Thermo Scientific Teflon/Rubber Laminated Discs (pkg. of 72)

Product #	Size	Fits These Containers
TS-12412	12 mm	0.3 and 1.0 ml Reacti-Vial Small Reaction Vials, 3.5 ml Screw Cap Septum Vials
TS-12418	18 mm	3 and 5 ml Reacti-Vial Small Reaction Vials
TS-12422	22 mm	25 and 40 ml Screw Cap Septum Vials

# Mininert Valves

Ideal for chemicals that deteriorate or evaporate in conventionally sealed containers.



They're easy to use. Push the green button to open, insert syringe needle and take sample, withdraw needle, then push red button to close. To change needle-seal septa, simply push the old septum out with a 1/8" diameter rod and push a new cylinder septum in. This is done with the valve closed to prevent exposure of contents.

Thermo Scientific Mininert Push-button Valves are highly dependable leak-tight closures for Screw Cap Septum Vials and other laboratory containers. Constructed of chemical-resistant Teflon coating, the valves provide an inert, high-pressure seal.

Mininert Valves are a superior replacement for rubber septum stoppers and ordinary screw caps. You can easily access the contents by inserting a syringe needle. A rubber gasket above the Teflon coating valve stem provides a seal for the needle when the valve is open. The seal prevents leakage and exposure of the contents during sampling.

# Mininert Valves are unique and practical seals for these containers:

- Screw Cap Septum Vials
- Thermo Scientific Hypo-Vial Sample Storage Vials
- Thermo Scientific Reacti-Vial Small Reaction Vials

Orderin	g Information			
Product #	Description	Fits these Containers	Size	Pkg. Size
TS-10135	Mininert Valves	3 and 5 ml Reacti-Vial Small Reaction Vials	20 mm	12/pkg
TS-10130	Mininert Valves	40 ml Screw Cap Septum Vials	27 mm	12/pkg

# Disc and Septa Compatibility Guide

Closure Type	Resealability	Recommended For Use With	Not Recommended For Use With
Teflon/Silicone Discs	Excellent	DMF, DMSO, organic solvents, pyridine, THF and silylation reactions	Strong corrosives, such as chlorosilanes
Teflon/Rubber Laminated Discs	Poor	Corrosives such as chlorosilanes, DMF, DMSO, organic solvents, pyridine and THF	Trifluoracetic anhydride

# Vacuum Hydrolysis Tubes

# For fast, effective protein and peptide hydrolysis

- The upper temperature limit of the Vacuum Hydrolysis Tubes is 260°C; however, do not heat the tubes greater than 100°C in an oven
- Vacuum Hydrolysis Tubes fit conveniently into Reacti-Block Aluminium Heating Blocks

Orderin	g Information			
Product #	Description	Fits these Containers	Size	Pkg. Size
TS-29570	Vacuum Hydrolysis	Tube	1mL	1 Each
TS-29571	Vacuum Hydrolysis	Tube	6mL	1 Each
TS-29572	Vacuum Hydrolysis	Tube	18mL	1 Each

Acetonitrile			14
Acetonitrile, HPLC Grade			42
Acylation Reagents		15-	18
Alkylation Reagents			
Aluminum Blocks, Reacti-Block			
Amino Acid Analysis			
Amino Acid Standard H			
N,O-Bis (trimethylsilyl)trifluoroacetam			
<i>N,O</i> -Bis (trimethylsilyl)trifluoroacetam			•
1% TMCS			8
10% TMCS			
BF <sub>2</sub> -Methanol			
Books, derivatization			
p-Bromophenacylate Reagent			
BSA			
BSTFA			
BSTFA 1% TMCS			
BSTFA + 10% TMCS			
Calibration Standard, HPLC			
Caps, Open-Top Screw			
Chromatography, HPLC			
Closures		20-	43
Discs	10	40	EΟ
Mininert Valves			
Open-Top Screw Caps			
Septa  Constant Boiling (6N) Hydrochloric Ac			
= -			
Dabsyl Chloride			
Derivatization Reagents, GC			
Derivatization Reagents, HPLC			
Detection Reagents, HPLC		.კე-	3/
4-Dimethylaminoazobenzene-			4.1
4-Sulfonyl Chloride			
Dimethylformamide		14,	42
N,N-Dimethylformamide dimethylace			
(Methylate Reagent)			
Dimethylsulfoxide			
Discs, Teflon Rubber	,		
Discs, Teflon/Silicone			
Dry Block Heaters			
Evaporators, Reacti-Vap			
FDAA, Marfey's Reagent			
Fluorinated Anhydrides			
Fluoracyl Imidazoles			17
1-Fluoro-2,4-dinitrophenyl-5-			
L-alanine amide			
Gas Chromatography			
Gas Chromatography Reagents		6-	25
Handbook of Analytical			
Derivatization Reactions			
Heating Blocks, Reacti-Block		46,	47

Heating Modules, Reacti-Therm	
Heptafluorobutyric Acid	26
Heptafluorobutyrylimidazole	17
Hexamethyldisilazane	11, 24
HFAA	
HFBI	17
HMDS	
HPLC	
Detection Reagents	29-31
HPLC Grade Solvents	
Ion Pair Reagents	
Hydrochloric Acid,	
Constant Boiling (6N)	39
Hydrolysis	
Hydrolysis Reagents	
Ion Pair Reagents	
Magnetic Stirrers	
Marfey's Reagent	
MBTFA	
MethElute Reagent	
Methylate Reagent	
N-Methyl-N-(tert-butyldimethylsilyl)	ZU
trifluoroacetamide	10
N-Methyl-bis (trifluoroacetamide)	
•	10
N-Methyl-	0
N-trimethylsilyltrifluoroacetamide	
Mininert Valves	
MOX Reagent	
MSTFA	
MSTFA + 1% TMCS	
MTBSTFA	
MTBSTFA + 1% TBDMCS	
Needles, Reacti-Vap Teflon Coated	
Ninhydrin	
Pentafluorobenzyl Bromide	
Pentafluoropropanol	
Pentafluoropropionic Acid Anhydride	
Peptide Retention Time Calibration Mix	
Peptide Retention Standard, HPLC	
Perfluoro acid anhydrides	
Perfluoroacylimidazoles	
PFAA	
PFBBr	
Phenylisothiocyanate	
PITC	
Protein Hydrolysis	38, 39
Pyridine	
Reacti-Block Aluminum Blocks	46-47
Reacti-Therm Heating Modules	44
Reacti-Therm Heating/ Stirring Modules	s44
Reacti-Vap Evaporators	44, 45

Reacti-Vap Teflon Coated Needles	
Reacti-Vial Magnetic Stirrers	
Reacti-Vial Small Reaction Vials	48
Reagents, Ion Pair	
Screw Cap Septum Vials	49
Screw Caps, Open-Top	48, 49
Septa	
Compatibility Guide	50
Septum Vials, Screw Cap	
Silicone Septa	
Silanes	
Siliconizing Fluids	
Silylation Grade Solvents	
Silylation Reagents	
Solvents (amino acid analysis)	
Solvents, HPLC	
Solvents, Silylation	
Standard, Peptide Retention	
Teflon	40
Needles	15
Rubber Laminated Discs	
Silicone Discs	
Valves	
Tetrahydrofuran	
TFAA	
TFAI	
Thermometer, Reacti-Therm	
TMCS	
TMPAH (MethElute Reagent)	
TMSI	
TNBSA	
Triethylamine	
Trifluoroacetic Acid	
Trifluoroacetic Acid Anhydride	
Trifluoroacetylimidazole	
Trimethylchlorosilane	
N-Trimethylsilylimidazole	
Trinitrobenzene Sulfonic Acid	35
Tri-Sil BP	13
Tri-Sil HTP	12
Tri-Sil Reagents	
Tri-Sil TBT	13
Tri-Sil TP	
Vacuum Hydrolysis Tubes	
Valves, Mininert	48, 49, 51
Vials	
Reacti-Vial Small Reaction	48
Screw Cap Septum Vials	49
	//2

TS-10130	48, 49, 51	TS-18821	44	TS-42799	23
TS-10135	48, 51	TS-18822	44	TS-42800	24
TS-12412	48, 49, 50	TS-18823	44	TS-42801	24
TS-12418	48, 50	TS-18824	44	TS-42855	24
TS-12422	48, 49, 50	TS-18825	44, 45	TS-44211	17
TS-12708	48, 49, 50	TS-18826	44, 45	TS-45950	12
TS-12712	48, 49, 50	TS-18827	45	TS-48882	17
TS-12713	49, 50	TS-18828	45	TS-48891	31
TS-12716	49, 50	TS-18829	45	TS-48895	31, 41
TS-12718	48, 50	TS-18830	45	TS-48910	9
TS-12720	50	TS-18831	45	TS-48911	9
TS-12722	48, 49, 50	TS-18832	45	TS-48913	9
TS-13019	49	TS-18833	45	TS-48914	9
TS-13028	49	TS-18914	46	TS-48915	9
TS-13043	49	TS-18915	46	TS-48920	10
TS-13080	49	TS-20062	14	TS-48927	10
TS-13097	48	TS-20088	39, 40	TS-48999	12
TS-13099	48	TS-20672	14	TS-49001	12
TS-13100	48	TS-20673	42	TS-49012	13
TS-13208	48, 49	TS-20684	14	TS-49016	13
TS-13215	48, 49	TS-20688	42	TS-49230	14
TS-13216	49	TS-21003	37	TS-49231	14
TS-13217	49	TS-21720	41	TS-49300	21
TS-13218	48, 49	TS-24308	39	TS-49301	21
TS-13219	48	TS-25003	26	TS-49350	20
TS-13220	48	TS-25104	42	TS-49370	19
TS-13221	48	TS-25108	27	TS-49700	16
TS-13222	48	TS-26922	40	TS-49701	16
TS-13223	48	TS-27530	14	TS-49703	16
TS-13225	48	TS-27860	14	TS-49704	16
TS-15012	25	TS-28901	28	TS-51101	42
TS-16000	47, 48	TS-28902	28	TS-51140	42
TS-16010	47, 48	TS-28903	28	TS-53101	27
TS-18782	45	TS-28904	28	TS-53104	26
TS-18784	45	TS-28905	27	TS-58220	20
TS-18786	45	TS-28997	35	TS-63162	18
TS-18801	46	TS-29570	51	TS-63163	18
TS-18802	46	TS-29571	51	TS-63164	18
TS-18803	46	TS-29572	51	TS-65191	18
TS-18804	46	TS-31700	43	TS-65192	18
TS-18806	47	TS-38828	8	TS-65193	18
TS-18807	47	TS-38829	8	TS-65195	15
TS-18809	47	TS-38830	8	TS-67363	18
TS-18810	47	TS-38831	8	TS-84769	11, 24
	46	TS-38832	8	TS-84770	11, 24
TS-18812	47	TS-38833	8	TS-88320	43
TS-18816	46	TS-38834	8	TS-88321	43
TS-18817	46	TS-38836	9	TS-88530	
TS-18818	46	TS-38838	9	TS-88623	10
TS-18819	46	TS-38839		TS-88625	
TS-18820	44	TS-38840	8	TS-88626	10

# Resources

# for Chromatographers

# Thermo Scientific Chromatography Columns and Consumables Catalog

This extensive catalog offers 590 pages of proven chromatography tools and product selection guides. Available online, with a robust search tool and optimized for your iPad®.

Visit www.thermoscientific.com/catalog



# Chromatography Resource Center

Our web-based resource center provides technical support, applications, technical tips and literature to help move your separations forward.

Visit www.thermoscientific.com/crc



# For more information visit: www.thermoscientific.com/crc

© 2012 Thermo Fisher Scientific Inc. All rights reserved. iPad is a registered trademark of Apple Inc., registered in the US and other countries. Dowex is a trademark of The Dow Chemical Company. Teflon is a registered trademark of E.I. du Pont de Nemours & Co., Inc. Cellosolve are registered trademarks of Union Carbide. Spheri-5 is a registered trademark of PerkinElmer. Mininert is a registered trademark of Dynatech Precision Sampling Corporation. Hycar is a registered trademark of Lubrizol Advance Materials. Luer-Lok is a registered trademark of Becton Dickinson and Company. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

USA and Canada +1 800 332 3331 France +33 (0)1 60 92 48 34 Germany +49 (0) 2423 9431 20 or 21 United Kingdom +44 (0)1928 534110 Japan +81 3 5826 1615 **China** +86 21 68654588 +86 10 84193588 +86 20 83145199 800 810 5118 **India** +91 22 6742 9494 +91 27 1766 2352 **Australia** 1 300 735 292 (free call domestic

India +91 22 6742 9494 +91 27 1766 2352 Australia 1 300 735 292 (free call domestic) New Zealand 0800 933 966 (free call domestic) All Other Enquiries +44 (0) 1928 534 050 **Technical Support** 

North America +1 800 332 3331 Outside North America +44 (0) 1928 534 440

