



Specification

Model	SH-DO-54FGE
Air Circulation	Forced
Temp Range	Ambient +10°C to 220°C (Ambient +18°F to 428°F)
Time Range	00.00 to 99HR 59MIN (MIN SEC) Selectable Digital Counter Timer
Temp & Time controller	Microprocess PID Membrane Touch Digital controller
Heater Capacity	800 W
Dimension Internal (WxDxH)	380×340×420(mm) 14.96×13.39×16.54(inch)
Dimension External (WxDxH)	531×612×760(mm) 20.91×24.09×29.92(inch)
Net Weight	44 Kg
Material Internal	Stainless Steel
Material External	Steel Plate With Powder Heating Coated
Shelves	2 Supplied (12 Maximum)
Electrical Requirements 120V, 50/60Hz, 1Φ	7.3 A
Capacity	54 L
Electrical Requirements 230V, 50/60Hz, 1Φ	3.5 A

Feature

- Automatic reset function incoming power cut
- Buzzer on when set time arrived
- Door open alarm
- Electric leakage breaker
- Over Temp Controller

Option

- Digital Recorder
- 2CH Hybrid Recorder