



Specification

Model	SH-FU-3MGE
Temp Controller	Digital PID Controller
Capacity	3 L
Heater Capacity	1100 W
Dimension Internal (WxDxH)	130×220×110(mm) 5.12×8.66×4.33(inch)
Dimension External (WxDxH)	300×378×405(mm) 11.81×14.88×15.94(inch)
Net Wight	19 Kg
Heater Element	KANTAL A-1 with CERAMIC FIBER
Insulation	Ceramic Board & Wool
Electrical Requirements 120V ,50/60Hz, 1Φ	10 A
Max Temp	1050°C (2192°F)
Electrical Requirements 230V ,50/60Hz, 1Φ	5A

Feature

Option
