



Specification

Model	SH-WB-6GAN
Demension Internal (WxDxH)	300×150×150(mm) 11.81×5.91×5.90(inch)
Demension External (WxDxH)	340×190×295(mm) 13.39×7.48×11.61(inch)
Capacity	6L
Temp Range	Ambient +5°C to 99°C (Ambient + 9°F to 210°F)
Time Range	00.00 to 99HR 59MIN (MIN SEC) Selectable Digital Display
Temp Controller	Temp Dial Controller
Net Weight	6 Kg
Material Internal	Stainless Steel
Material External	Steel Plate With Powder Heating Coated
Material Lid	Stainless Steel
Heater Capacity	600 W
Circulation Pump	None
Drain Valve	None
Electrical Requirements 120V ,50/60Hz, 1Φ	5.1 A
Electrical Requirements 230V,50/60Hz, 1Φ	2.7 A

Feature

- Automatic reset function incoming power cut
- Buzzer on when set time arrived
- Electric leakage circuit breaker

Option

- Heater cover
- Bath Lid