



Specification

Model	SH-CH-150U1
Dimension Internal (WxDxH)	500x500x600(mm) 19.69x19.69x23.62(inch)
Dimension External (WxDxH)	920x995x1555(mm) 36.22x39.17x61.22(inch)
Capacity	150 L
Heater Capacity	2 KW
Temp Uniformity	±0.5 °C
Temp Range	-25 °C ~ 120 °C (-13 °F ~ 248 °F)
Cool Down Time	23 °C ~ -20 °C (73.4 °F ~ -4 °F) About 40min
Heating up Time	-20 °C ~ 100 °C (-4 °F ~ 212 °F) About 40min
Temp Controller	Microprocessor Programmable Controller Touch Screen Type
Material Internal	Stainless Steel
Material External	Steel Plate with Powder Heating Coated
Electrical Requirements 230V ,50/60Hz, 1Φ	14.5 A
Shelves	Movable stainless wire shelves 2 Supplied (17 Maximum)
Temp Accuracy & Stability	±0.1 °C
Electrical Requirements 230V,50/60Hz, 3Φ	N/A
Weight	250 Kg

Feature

- Automatic reset function incoming power cut
- Over temp protection
- Electric leakage breaker
- Door open lamp display
- Fan motor over load protection
- Refrigerator over load protection

Option

- Digital recorder
- 2CH Hybrid recorder
- Max temperature 150 °C