



Microplate Washers

Expand functionality beyond the standard unit

Our reliable and easy to use Thermo Scientific™ Wellwash™ microplate strip washers provide secure washing performance for routine and research ELISA applications. Wellwash washers feature a large color display with graphical user interface, local language versions, and non-pressurized bottles for optimal ease-of-use, convenience, and safety.

Which Wellwash product is best for your needs?



[Wellwash](#)



[Wellwash Versa](#)

Wash head	1x8, 1x12	two 1x8, two 1x12, 2x8 cell wash head
Wash bottles	1 bottle (2 liter)	3 bottles (two 4 liter and one 2 liter)
Plate formats	96	96, 384
Cell washing	No	Yes
Robot compatibility	No	Yes

[Microplate reader brochure>](#)

